

Andrea Brazdova, Oumsaad Naas, Nicolas Visez, Jean-Pierre Sutra, H el ene S en echal et Pascal Poncet

H opital d'Enfants Armand Trousseau, Laboratoire de Biochimie, Equipe «Allergie & Environnement», 26 avenue du Dr Arnold Netter, 75012 - Paris.
Universit e de Lille, Laboratoire de Physico-chimie des

Processus de Combustion et de l'Atmosph ere, Lille. Institut Pasteur, d epartement Infection & Epid emiologie, Paris. 2015

Overview

- **Keywords:** Allergen, IgE immunoreactivity, pollen homogenization, hypersensitivity community, DNA extraction
- **Aim of the study:** Identification of a fast method for protein extraction from pollen grains
- **Application:** Western blot analysis
- **Sample name:** Birch, Nettle, Wall Pellitory pollens
- **Sample type:** Pollen
- **Material:** FastPrep-24™ 5G instrument, CoolPrep adapter, 2 mL Lysing Matrix C & E tubes
- **Buffer:** PBS

Protocol and Parameters

Incubation Method

- 1. Add 50 mg of pollen and 500 µL of PBS in a tube
- 2. Place the tube in a shaker for 18 hours in cold room
- 3. Centrifuge the suspension 20 mins at 18,000 x g, 4°C
- 4. Keep the supernatant at -20°C prior to analysis

Grinding Method

- 1. Add 50 mg of pollen and 500 µL of PBS to a 2 mL Lysing Matrix C or E tube.
- 2. Load Lysing Matrix tubes in a CoolPrep Adapter containing dry ice.
- 3. Process with the FastPrep-24 5G: 40 sec at a speed setting of 6.0 m/s.
- 4. Centrifuge the Lysing Matrix tubes 20 mins at 18,000 x g, 4°C to pellet debris.
- 5. Keep the supernatant at -20°C prior to analysis

CASE STUDY

Conclusion


Protein extraction from pollen samples with the FastPrep-24™ 5G showed to be highly effective compared to the standard method based on overnight incubation. The effectiveness of the FastPrep method is quantitative, resulting in higher protein yields, as well as qualitative, as evidenced by a wide variety composition of protein extracts. The FastPrep system is a powerful tool to rapidly generate protein extracts with high reproducibility, ready for electrophoresis (SDS-PAGE) analysis. IgE immunoreactivity is conserved in protein extracted with the FastPrep-24™ 5G instrument.

Total destruction of the pollen grain structure with FastPrep-24™ 5G instrument and Lysing Matrix C



Optical microscope observation of pine pollen (x200) before (left) and after grinding (right) with the FastPrep-24™ 5G System.

Up to 2 mg/mL of protein extracted with FastPrep-24™ 5G System



Comparison of 8 pine pollen protein extracts obtained by standard or FastPrep method. Experiment were repeated 4 times using 4 different pollen batches. Protein concentration was determined using Bradford assay.

■ FastPrep-24™
■ Incubation/Rotation

MP Biomedicals

Americas: 800.854.0530 | custserv@mpbio.com

Europe: 00800.7777.9999 | custserv.eur@mpbio.com

APAC: 65.6775.0008 | enquiry_ap@mpbio.com

LEARN MORE!

www.mpbio.com