



Taq Core Kit 10

Cat. # 11EPTQK10B

Taq Core Kits 10 contain all reagents required for PCR but as separated items:

- **Taq DNA Polymerase**, a highly purified thermostable recombinant polymerase, with a 5'-3' exonuclease activity, but no 3'-5' exonuclease activity. *Taq* DNA Polymerase resists to prolonged incubations at 95°C.
- Highly purified dNTPs in a ready-to-use mastermix at 10mM each, guaranteed for PCR.
- Optimised standard PCR buffers 10xC provided with and without MgCl₂, highly efficient in most PCR conditions.
- MgCl₂ solution at 25mM as a separate item for peculiar optimisation of MgCl₂ concentration.

Reaction conditions:

Composition of PCR buffer with MgCl₂ diluted to 1xC: 10 mM Tris-HCl pH 9.0 (25°C); 50 mM KCl; 1.5 mM MgCl₂; 0.1% Triton X100; 0.2 mg/ml BSA.

Taq DNA polymerase is supplied in a formulation containing 50% glycerol.

Composition of storage buffer: 20mM Tris-HCl pH 8.0; 100mM KCl; 0.1mM EDTA; 1mM Dithiothreitol; 0.5 % Tween 20; 0.5% Nonidet P40; 50% glycerol.

The following information is provided as general advice. When setting up multiple reactions, the preparation of a reaction mastermix containing all non variable components is recommended. That will reduce risks of cross contamination and errors when pipetting low volumes.

Contents	
No. of Tubes	Components
5	1000U <i>Taq</i> DNA Polymerase at 5 U/μl
14	dNTPs mix at 10mM each
14	PCR buffer 10xC with MgCl ₂
5	PCR buffer 10xC without MgCl ₂
3	MgCl ₂ at 25 mM

Store at -20°C.

For research use only.

Reagent	20 μl Reaction	50 μl Reaction	Final Concentration
Sterilized water	Add to 20μl	Add to 50μl	
PCR buffer 10xC	2 μl	5 μl	1xC
dNTPs mix 10 mM each	0.4 μl	1 μl	200 μM each
Primer A ***	x μl	x μl	0.2 to 0.5 μM
Primer B ***	x μl	x μl	0.2 to 0.5 μM
DNA Template *	x μl	x μl	
<i>Taq</i> DNA Pol 5 U/μl	0.2 μl to 0.3 μl	0.2 μl to 0.3 μl	1 to 1.5 U/rxn **

* Suggested amount of DNA template: 100 pg of plasmid or phage DNA; 10 to 100 ng of genomic DNA; 0.2 to 30 ng of cDNA.

** If PCR inhibitors remain in the reaction mix, higher amounts of *Taq* DNA Pol may be necessary (2-3 U/rxn) to ensure optimal PCR results. Same advice for amplification of fragments up to 3 – 4 kb.

*** Ideally, prepare a 10x solution containing both PCR primers.

Optimisation of Mg concentration:

Most PCR applications are optimal with 1.5mM MgCl₂. If optimisation of MgCl₂ concentration is necessary, use the PCR buffer without MgCl₂ and the separate MgCl₂ vial at 25mM.

The recommended range of MgCl₂ concentration is between 1mM and 4mM.

Final concentration of MgCl ₂ (mM)	1	1.25	1.5	1.75	2	2.5	3	4
Volume (μl) of 25mM MgCl ₂ in a 20 μl PCR reaction	0.16	1	1.2	1.4	1.6	2	2.4	3.2
Volume (μl) of 25mM MgCl ₂ in a 50 μl PCR reaction	2	2.5	3	3.5	4	5	6	8

Guidelines for PCR program:

Cycle step	Temperature	Time	Cycles
Initial denaturation (1)	93-95°C	2-4 min	1
Denaturation	93 °C	30 - 60 sec	20 -40
Annealing (2)	(T _m – 5°C)	30 - 60 sec	
Elongation (3)	72°C	1 min /kb	
Final elongation (4)	72°C	10 min	1
End of PCR assay	+ 4°C / -20°C	Hold / Store	

(1) For GC rich templates, extend to 10 min, but add *Taq* DNA Pol after initial denaturation or keep maximum 5 min at 95°C by adding up to 10% DMSO to reaction mix.

(2) The annealing temperature depends on the length and the T_m (melting temperature) of the primers. It should be placed 5°C below the lowest T_m of the pair of primers. Avoid annealing temperatures under 48°C.

(3) Elongation time will depend on the length of the amplified template.

For amplifications less than 1kb, 15 to 30 sec could be enough.

(4) A final elongation step is recommended to fill-in the incompletely amplified fragments and to add the A residues to the 3'ends of PCR products, (useful for T/A cloning protocols).

Characteristics of dNTPs:

- dNTPs Mix 10 mM is a premixed solution in a TE Buffer (Tris-HCl 10mM, EDTA 1mM, pH7,5), containing the sodium salts of 4 nucleotides (dATP, dCTP, dGTP, dTTP), each at a concentration of 10mM.
- The individual dNTPs, PCR grade, are guaranteed of high purity as measured by HPLC (dATP, dCTP, dGTP, and dTTP HPLC areas % >99%).
- The solution is specially manufactured and tested for application in PCR, RT-PCR, qPCR, DNA Labeling and sequencing.

Guidelines for primer sequence design:

- T_m (melting temperature) of both primers should be equivalent and not exceed 70°C.
- Annealing temperature (T_a) should be about 5°C below the lowest T_m of the pair of primers.
- One primer length will depend upon its base content and the T_m of its partner.
- The likelihood of annealing to sequences other than the chosen target has to be very low. The size of the primers must be over 16-18 bases for human or animal templates and over 20-22 bases for plants. The length should not exceed 28-30 bases.
- Final base composition should always be 50-60% (G+C).
- 3' ends of primers should not be complementary to avoid primer dimers.
- Avoid 3 G or C residues in the 3' end region of the primer, as it may result in non specific annealing.
- Ability to form secondary structures such as hairpins should be avoided.

Quality control:

Unit assay conditions:

One unit of DNA polymerase is the amount of enzyme required to catalyse the incorporation of 10 nanomoles of nucleosides into a DE81 absorbable product within 30 min at 74°C under assay conditions.

Absence of contamination: No nickases, endonucleases, 3' exonuclease and ribonucleases activities detected.

PCR assay on human genomic DNA: *Taq* DNA Polymerase, PCR buffers and dNTPs mix, are analysed by PCR to check batch to batch equivalence. Specific PCR products of 400 bp using a human β -globin gene as a template and decreasing amounts of both DNA template and *Taq* DNA polymerase are obtained in all assays.

Troubleshooting:

1. Little or no amplification observed:

- Increase the elongation time on the basis of 60 sec/kb.
- Increase the amount of *Taq* DNA Polymerase directly up to 2U/rxn or in 0.5 U steps.
- Increase the concentration of MgCl₂ in 0.25mM steps (see the table for optimisation of Mg concentration).
- Check the T_m of the primers, their quality and their length (>18 bases for animal material and >20-22 bases for plants).
- Avoid primer secondary structures, like dimers, hairpins.
- Avoid any mismatch or point mutations in the 3' region of the primers.
- Avoid 3G or 3C residues in the 3' end region of the primer.
- Both primers might be present at the same final concentration.
- Prevent primers degradation by making aliquots and storing them at -20°C.
- Lower the annealing temperature by 2°C to 5°C steps.
- If annealing temperature (T_a) is close to the elongation temperature, perform PCR by decreasing "T_a" every cycle.
- Increase the number of cycles by 5 (if < 35).
- Make sure that the final elongation of 10 min at 72°C has been programmed.
- If templates are GC rich or contain difficult secondary structures, add DMSO up to 10% by 2% steps. Reduce enzyme concentration (0.5 U/rxn).
- Increase the quantity of the template, and/or the primers.
- Verify the quality of your DNA template and check for possible degradation by loading onto an agarose gel.
- When using cDNA as a template, check the quality of your RNA and consider the performance of RT step. Load an aliquot of cDNA onto an agarose gel.
- The volume of RT-reaction (cDNA) should not exceed 10% of the PCR reaction volume.
- Consider the presence of inhibitory factors.
- Use thin-wall PCR tubes.

2. Multiple bands or smearing observed:

- Decrease the amount of *Taq* DNA Polymerase in 0.5 U steps.
- Check the T_m of the primers, their length, their quality.
- Check the homologous sequences of the primers to the template.
- Verify the absence of secondary structures, like hairpins.
- Both primers might be present at the same final concentration.
- Increase the annealing temperature by 2°C to 5°C steps.
- Lengthen the homologous sequence of the primer to the template.
- Reduce the number of cycles by 5.
- Reduce the concentration of DNA template and/or primers.
- Use serial dilution of DNA template.

For further Information contact us:

MP Biomedicals Europe

Tel: 00800.7777.9999

Fax: 00800.6666.8888

custserv.eur@mpbio.com

MP Biomedicals North America

Tel: 800 854 05 30

Fax: 800 934 69 99

custserv@mpbio.com

www.mpbio.com