



## Taq Core Kits 25

Cat. # 11EPTQK251

**Taq Core Kits 25** contain all reagents required for PCR but as separated items:

- *Taq* DNA Polymerase, a highly purified thermostable recombinant polymerase, with a 5'-3' exonuclease activity, but no 3'-5' exonuclease activity. *Taq* DNA Polymerase resists to prolonged incubations at 95°C.
- Highly purified dNTPs in a ready-to-use mastermix at 25 mM each, guaranteed for PCR.
- Optimised standard PCR buffers 10xC provided with and without MgCl<sub>2</sub>, highly efficient in most PCR conditions.
- MgCl<sub>2</sub> solution at 25mM as a separate item for peculiar optimisation of MgCl<sub>2</sub> concentration.

### Reaction conditions:

Composition of PCR buffer with MgCl<sub>2</sub> diluted to 1xC: 10 mM Tris-HCl pH 9.0 (25°C); 50 mM KCl; 1.5 mM MgCl<sub>2</sub>; 0.1% Triton X100; 0.2 mg/ml BSA.

*Taq* DNA polymerase is supplied in a formulation containing 50% glycerol.

Composition of storage buffer: 20mM Tris-HCl pH 8.0; 100mM KCl; 0.1mM EDTA; 1mM Dithiothreitol; 0.5% Tween 20; 0.5% Nonidet P40; 50% glycerol.

The following information is provided as general advice. When setting up multiple reactions, the preparation of a reaction mastermix containing all non variable components is recommended. That will reduce risks of cross contamination and errors when pipetting low volumes.

Contents	
No. of Tubes	Components
1	250U <i>Taq</i> DNA Polymerase at 5 U/μl
1	dNTPs mix at 25mM each
2	PCR buffer 10xC with MgCl <sub>2</sub>
1	PCR buffer 10xC without MgCl <sub>2</sub>
1	MgCl <sub>2</sub> at 25 mM

Reagent	20 μl Reaction	50 μl Reaction	Final Concentration
Sterilized water	Add to 20μl	Add to 50μl	
PCR buffer 10xC	2 μl	5 μl	1xC
dNTPs mix 25 mM each	0.16 μl	0.4 μl	200 μM each
Primer A ***	x μl	x μl	0.2 to 0.5 μM
Primer B ***	x μl	x μl	0.2 to 0.5 μM
DNA Template *	x μl	x μl	
<i>Taq</i> DNA Pol 5 U/μl	0.2 μl to 0.3 μl	0.2 μl to 0.3 μl	1 to 1.5 U/rxn **

\* Suggested amount of DNA template: 100 pg of plasmid or phage DNA; 10 to 100 ng of genomic DNA; 0.2 to 30 ng of cDNA.

\*\* If PCR inhibitors remain in the reaction mix, higher amounts of *Taq* DNA Pol may be necessary (2-3 U/rxn) to ensure optimal PCR results. Same advice for amplification of fragments up to 3 – 4 kb.

\*\*\*Ideally, prepare a 10x solution containing both PCR primers.

### Optimisation of Mg concentration:

Most PCR applications are optimal with 1.5mM MgCl<sub>2</sub>. If optimisation of MgCl<sub>2</sub> concentration is necessary, use the PCR buffer without MgCl<sub>2</sub> and the separate MgCl<sub>2</sub> vial at 25mM.

The recommended range of MgCl<sub>2</sub> concentration is between 1mM and 4mM.

Final concentration of MgCl <sub>2</sub> (mM)	1	1.25	1.5	1.75	2	2.5	3	4
Volume (μl) of 25mM MgCl <sub>2</sub> in a 20 μl PCR reaction	0.16	1	1.2	1.4	1.6	2	2.4	3.2
Volume (μl) of 25mM MgCl <sub>2</sub> in a 50 μl PCR reaction	2	2.5	3	3.5	4	5	6	8

### Guidelines for PCR program:

Cycle step	Temperature	Time	Cycles
Initial denaturation (1)	93-95°C	2-4 min	1
Denaturation	93 °C	30 - 60 sec	20 - 40
Annealing (2)	(T <sub>m</sub> – 5°C)	30 - 60 sec	
Elongation (3)	72°C	1 min /kb	
Final elongation (4)	72°C	10 min	1
End of PCR assay	+ 4°C / -20°C	Hold / Store	

(1) For GC rich templates, extend to 10 min, but add *Taq* DNA Pol after initial denaturation or keep maximum 5 min at 95°C by adding up to 10% DMSO to reaction mix.

(2) The annealing temperature depends on the length and the T<sub>m</sub> (melting temperature) of the primers. It should be placed 5°C below the lowest T<sub>m</sub> of the pair of primers. Avoid annealing temperatures under 48°C.

(3) Elongation time will depend on the length of the amplified template.

For amplifications less than 1kb, 15 to 30 sec could be enough.

(4) A final elongation step is recommended to fill-in the incompletely amplified fragments and to add the A residues to the 3'ends of PCR products, (useful for T/A cloning protocols).

Store at –20°C.

For research use only.

### Characteristics of dNTPs:

- dNTPs Mix 25mM is a premixed solution in a TE Buffer (Tris-HCl 10mM, EDTA 1mM, pH7.5), containing the sodium salts of 4 nucleotides (dATP, dCTP, dGTP, dTTP), each at a concentration of 25mM.
- The individual dNTPs, PCR grade, are guaranteed of high purity as measured by HPLC (dATP, dCTP, dGTP, and dTTP HPLC areas % >99%).
- The solution is specially manufactured and tested for application in PCR, RT-PCR, qPCR, DNA Labelling and sequencing.

### Guidelines for primer sequence design:

- T<sub>m</sub> (melting temperature) of both primers should be equivalent and not exceed 70°C.
- Annealing temperature (T<sub>a</sub>) should be about 5°C below the lowest T<sub>m</sub> of the pair of primers.
- One primer length will depend upon its base content and the T<sub>m</sub> of its partner.
- The likelihood of annealing to sequences other than the chosen target has to be very low. The size of the primers must be over 16-18 bases for human or animal templates and over 20-22 bases for plants. The length should not exceed 28-30 bases.
- Final base composition should always be 50-60% (G+C).
- 3' ends of primers should not be complementary to avoid primer dimers.
- Avoid 3 G or C residues in the 3' end region of the primer, as it may result in non specific annealing.
- Ability to form secondary structures such as hairpins should be avoided.

### Quality control:

#### Unit assay conditions:

One unit of DNA polymerase is the amount of enzyme required to catalyse the incorporation of 10 nanomoles of nucleosides into a DE81 absorbable product within 30 min at 74°C under assay conditions.

Absence of contamination: No nickase, endonuclease, 3' exonuclease, ribonuclease activities detected.

PCR assay on human genomic DNA: *Taq* DNA Polymerase, PCR buffers and dNTPs mix, are analysed by PCR to check batch to batch uniformity. Specific PCR products of 400 bp using a human β-globin gene as a template and decreasing amounts of both DNA template and *Taq* DNA polymerase are obtained in all assays.

### Troubleshooting:

#### 1. Little or no amplification observed:

- Increase the elongation time on the basis of 60 sec/kb.
- Increase the amount of *Taq* DNA Polymerase directly up to 2U/rxn or in 0.5 U steps.
- Increase the concentration of MgCl<sub>2</sub> in 0.25mM steps (see the table for optimisation of Mg concentration).
- Check the T<sub>m</sub> of the primers, their quality and their length (>18 bases for animal material and >20-22 bases for plants).
- Avoid secondary structure primers, like dimers, hairpins.
- Avoid any mismatch or point mutations in the 3' region of the primers.
- Avoid 3G or 3C residues in the 3' end region of the primer.
- Both primers might be present at the same final concentration.
- Prevent primers degradation by making aliquots and storing them at -20°C.
- Lower the annealing temperature by 2°C to 5°C steps.
- If annealing temperature (T<sub>a</sub>) is close to the elongation temperature, perform PCR by decreasing "T<sub>a</sub>" every cycle.
- Increase the number of cycles by 5 (if <35).
- Make sure that the final elongation of 10 min at 72°C has been programmed.
- If templates are GC rich or contain difficult secondary structures, add DMSO up to 10% by 2% steps. Reduce enzyme concentration (0.5 U/rxn).
- Increase the quantity of the template, and/or the primers.
- Verify the quality of your DNA template and check for possible degradation onto an agarose gel.
- When using cDNA as a template, check the quality of your RNA and consider the performance of RT step. Load an aliquot of cDNA onto an agarose gel.
- The volume of RT-reaction (cDNA) should not exceed 10% of the PCR reaction volume.
- Consider the presence of inhibitory factors.
- Use thin-wall PCR tubes.

#### 2. Multiple bands or smearing observed:

- Decrease the amount of *Taq* DNA Polymerase in 0.5 U steps.
- Check the T<sub>m</sub> of the primers, their length and their quality.
- Check the homologous sequences of the primers to the template.
- Verify the absence of secondary structures, like hairpins.
- Both primers might be present at the same final concentration.
- Increase the annealing temperature by 2°C to 5°C steps.
- Lengthen the homologous sequence of the primer to the template.
- Reduce the number of cycles by 5.
- Reduce the concentration of DNA template and/or primers.
- Use serial dilution of DNA template.

### Other PCR buffers available from MP Biomedicals:

- **PCR buffer XD** formulated without BSA and without detergents, to work on plant or animal material or using further dHPLC.
- **(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> PCR buffer** containing ammonium sulphate to improve efficiency and specificity of some applications.
- **Direct loading buffer** with a densifying agent and a red purple dye that allows direct monitoring on agarose gels.

For further Information contact us:

#### **MP Biomedicals Europe**

Tel: 00800.7777.9999

Fax: 00800.6666.8888

[custserv.eur@mpbio.com](mailto:custserv.eur@mpbio.com)

#### **MP Biomedicals North America**

Tel: 800 854 05 30

Fax: 800 934 69 99

[custserv@mpbio.com](mailto:custserv@mpbio.com)

[www.mpbio.com](http://www.mpbio.com)