



Taq Core Kit 10

Cat. #
11EPTQK300

Pack size:
1000 U

| Contents | |
|--------------|---|
| No. of Tubes | Components |
| 1 | 1000U <i>Taq</i> DNA Polymerase at 5 U/ μ l |
| 3 | dNTPs mix at 10 mM each |
| 4 | PCR buffer 10xC with MgCl ₂ |
| 3 | PCR buffer 10xC without MgCl ₂ |
| 2 | MgCl ₂ at 25 mM |

Store at **-20°C**.

For research use only.

Taq Core Kits 10 contain all reagents required for PCR but as separated items:

- *Taq* DNA Polymerase, a highly purified thermostable recombinant polymerase, with a 5'-3' exonuclease activity, but no 3'-5' exonuclease activity. *Taq* DNA Polymerase resists to prolonged incubations at 95°C.
- Highly purified dNTPs in a ready-to-use mastermix at 10 mM each, guaranteed for PCR.
- Optimised standard PCR buffers 10xC provided with and without MgCl₂, highly efficient in most PCR conditions.
- MgCl₂ solution at 25 mM as a separate item for peculiar optimisation of MgCl₂ concentration

Reaction conditions:

Composition of PCR buffer with MgCl₂ diluted to 1xC: 10 mM Tris-HCl pH 9.0 (25°C); 50 mM KCl; 1.5 mM MgCl₂; 0.1% Triton X100; 0.2 mg/ml BSA.

Taq DNA polymerase is supplied in a formulation containing 50% glycerol.

Composition of storage buffer: 20 mM Tris-HCl pH 8.0; 100 mM KCl; 0.1 mM EDTA; 1 mM Dithiothreitol; 0.5 % Tween 20; 0.5 % Nonidet P40; 50 % glycerol.

The following information is provided as general advice. When setting up multiple reactions, the preparation of a reaction mastermix containing all non variable components is recommended. That will reduce risks of cross contamination and errors when pipetting low volumes.

| Reagent | 20 μ l Reaction | 50 μ l Reaction | Final Concentration |
|---------------------------------|----------------------------|----------------------------|---------------------|
| Sterilized water | Add to 20 μ l | Add to 50 μ l | |
| PCR buffer 10xC | 2 μ l | 5 μ l | 1xC |
| dNTPs mix 10 mM each | 0.4 μ l | 1 μ l | 200 μ M each |
| Primer A *** | x μ l | x μ l | 0.2 to 0.5 μ M |
| Primer B *** | x μ l | x μ l | 0.2 to 0.5 μ M |
| DNA Template * | x μ l | x μ l | |
| <i>Taq</i> DNA Pol 5 U/ μ l | 0.2 μ l to 0.3 μ l | 0.2 μ l to 0.3 μ l | 1 to 1.5 U/rxn ** |

* Suggested amount of DNA template: 100 pg of plasmid or phage DNA; 10 to 100 ng of genomic DNA; 0.2 to 30 ng of cDNA.

** If PCR inhibitors remain in the reaction mix, higher amounts of *Taq* DNA Pol may be necessary (2-3 U/rxn) to ensure optimal PCR results. Same advice for amplification of fragments up to 3 – 4 kb.

*** Ideally, prepare a 10x solution containing both PCR primers.

Optimisation of Mg concentration:

Most PCR applications are optimal with 1.5 mM MgCl₂. If optimisation of MgCl₂ concentration is necessary, use the PCR buffer without MgCl₂ and the separate MgCl₂ vial at 25 mM. The recommended range of MgCl₂ concentration is between 1 and 4 mM.

| | | | | | | | | |
|---|------|------|-----|------|-----|-----|-----|-----|
| Final concentration of MgCl ₂ (mM) | 1 | 1.25 | 1.5 | 1.75 | 2 | 2.5 | 3 | 4 |
| Volume (μ l) of 25 mM MgCl ₂ in a 20 μ l PCR reaction | 0.16 | 1 | 1.2 | 1.4 | 1.6 | 2 | 2.4 | 3.2 |
| Volume (μ l) of 25 mM MgCl ₂ in a 50 μ l PCR reaction | 2 | 2.5 | 3 | 3.5 | 4 | 5 | 6 | 8 |

Guidelines for PCR program:

| Cycle step | Temperature | Time | Cycles |
|--------------------------|------------------------|--------------|--------|
| Initial denaturation (1) | 93-95°C | 2-4 min | 1 |
| Denaturation | 93 °C | 30 - 60 sec | 20 -40 |
| Annealing (2) | (T _m – 5°C) | 30 - 60 sec | |
| Elongation (3) | 72°C | 1 min /kb | |
| Final elongation (4) | 72°C | 10 min | 1 |
| End of PCR assay | + 4°C / -20°C | Hold / Store | |

(1) For GC rich templates, extend to 10 min, but add *Taq* DNA Pol after initial denaturation or keep maximum 5 min at 95°C by adding up to 10% DMSO to reaction mix.

(2) The annealing temperature depends on the length and the T_m (melting temperature) of the primers. It should be placed 5°C below the lowest T_m of the pair of primers. Avoid annealing temperatures under 48°C.

(3) Elongation time will depend on the length of the amplified template.

For amplifications, less than 1kb, 15 to 30 sec could be enough.

(4) A final elongation step is recommended to fill-in the incompletely amplified fragments and to add the A residues to the 3'ends of PCR products, (useful for T/A cloning protocols).

Characteristics of dNTPs:

- dNTPs Mix 10 mM (Cat# NTPMX100) is a premixed solution in a TE Buffer (Tris-HCl 10mM, EDTA 1mM, pH7.5), containing the sodium salts of 4 nucleotides (dATP, dCTP, dGTP, dTTP), each at a concentration of 10 mM
- The individual dNTPs, PCR grade, are guaranteed of high purity as measured by HPLC (dATP, dCTP, dGTP, and dTTP HPLC areas % > 99%).
- The solution is specially manufactured and tested for application in PCR, RT-PCR, qPCR, DNA Labeling and sequencing.

Guidelines for primer sequence design:

- T_m (melting temperature) of both primers should be equivalent and not exceed 70°C
- Annealing temperature (T_a) should be about 5°C below the lowest T_m of the pair of primers.
- One primer length will depend upon its base content and the T_m of its partner.
- The likelihood of annealing to sequences other than the chosen target has to be very low. The size of the primers must be over 16-18 bases for human or animal templates and over 20-22 bases for plants. The length should not exceed 28-30 bases.
- Final base composition should always be 50-60% (G+C)
- 3' ends of primers should not be complementary to avoid primer dimers
- Avoid 3 G or C residues in the 3' end region of the primer, as it may result in non-specific annealing
- Ability to form secondary structures such as hairpins should be avoided.

Quality control:

Unit assay conditions:

One unit of DNA polymerase is the amount of enzyme required to catalyse the incorporation of 10 nanomoles of nucleosides into a DE81 adsorbable product within 30 min at 74°C under assay conditions.

Absence of contamination: No nickase, endonuclease, 3' exonuclease, ribonuclease activities detected.

PCR assay on human genomic DNA: *Taq* DNA Polymerase, PCR buffers and dNTPs mix, are analysed by PCR to check batch to batch equivalence. Specific PCR products of 400 bp using a human β -globin gene as a template and decreasing amounts of both DNA template and *Taq* DNA polymerase are obtained in all assays.

Troubleshooting:

1. Little or no amplification observed:

- Increase the elongation time on the basis of 60 sec/kb
- Increase the amount of *Taq* DNA Polymerase directly up to 2U/rxn or in 0.5 U steps
- Increase the concentration of MgCl₂ in 0.25 mM steps (see the table for optimisation of Mg concentration)
- Check the T_m of the primers, their quality and their length (> 18 bases for animal material and > 20-22 bases for plants)
- Avoid primer secondary structures, like dimers, hairpins
- Avoid any mismatch or point mutations in the 3' region of the primers
- Avoid 3 G or 3C residues in the 3' end region of the primer
- Both primers might be present at the same final concentration
- Prevent primers degradation by making aliquots and storing them at -20°C
- Lower the annealing temperature by 2°C to 5°C steps
- If annealing temperature (T_a) is close to the elongation temperature, perform PCR by decreasing "T_a" every cycle.
- Increase the number of cycles by 5 (if < 35)
- Make sure that the final elongation of 10 min at 72°C has been programmed
- If templates are GC rich or contain difficult secondary structures, add DMSO up to 10% by 2% steps. Reduce enzyme concentration (0.5 U/rxn).
- Increase the quantity of the template, and/or the primers
- Verify the quality of your DNA template and check for possible degradation by loading onto an agarose gel
- When using cDNA as a template, check the quality of your RNA and consider the performance of RT step. Load an aliquot of cDNA onto an agarose gel
- The volume of RT-reaction (cDNA) should not exceed 10% of the PCR reaction volume
- Consider the presence of inhibitory factors
- Use thin-wall PCR tubes

2. Multiple bands or smearing observed:

- Decrease the amount of *Taq* DNA Polymerase in 0.5 U steps
- Check the T_m of the primers, their length, their quality
- Check the homologous sequences of the primers to the template
- Verify the absence of secondary structures, like hairpins
- Both primers might be present at the same final concentration
- Increase the annealing temperature by 2°C to 5°C steps
- Lengthen the homologous sequence of the primer to the template
- Reduce the number of cycles by 5
- Reduce the concentration of DNA template and/or primers
- Use serial dilution of DNA template

Other PCR buffers available from MP Biomedicals:

- **PCR buffer XD** formulated without BSA and without detergents, to work on vegetal or beef material or using further dHPLC
- **(NH₂)SO₄ PCR buffer** containing ammonium sulphate to improve efficiency and specificity of some applications
- **Direct loading buffer** with a densifying agent and a red purple dye allowing direct loading after cycling without new handle.

For further Information contact us:

MP Biomedicals Europe

Tel: 00800.7777.9999

Fax: 00800.6666.8888

custserv.eur@mpbio.com

MP Biomedicals North America

Tel: 800 854 05 30

Fax: 800 934 69 99

custserv@mpbio.com

www.mpbio.com