

MagBeads FastDNA[®] Kit for Feces

Ready-to-Use for MPure-32[™] aNAP System



Cat. No.: 117033200 (96 Preps)

Quick-Start Manual

Revision 1.0 June 2022

Notes before starting

- ❑ Expect precipitation in Pre-Wash Buffer and Lysis Buffer F1, warm the solutions to 55 °C will dissolve the precipitate.
- ❑ Vortex the sample in a Lysing Matrix E tube at the maximum speed for 10 mins if a FastPrep[®] Instrument is unavailable. Secure samples on the vortex through an adapter to ensure homogenization.
- ❑ Centrifugation speed stated in the manual will be a guideline, use the maximum speed available if 14,000 x g is not feasible.

Automation Extraction

1. Add 30-300 mg fecal sample to **Lysing Matrix E** tube and resuspend it with 1000 µL of **Pre-Wash Buffer**. Vortex for 10-15 seconds.

Note: Ensure there is at least 200 µL of empty space in the tube after adding the sample.

2. Centrifuge samples at 14,000 rpm for 3-5 mins and remove supernatant.

3. Add 490 µL **Lysis Buffer F1**, 60 µL **Lysis Buffer F2** and 10 µL **RNase A Solution** to sample in Lysing Matrix E tube. Shake the tube to mix well.

4. Homogenize sample in FastPrep[®] Instrument for 40 seconds at a speed of 5.0 m/s.

Note: The speed and time can be changed according to different fecal samples. Vortex 5-10 mins at maximum speed if FastPrep[®] Instrument is not available.

5. Centrifuge for 5 mins at 14,000 x g to pellet down the debris at the bottom of tube.

6. Transfer 500 µL supernatant to a clean 2.0 mL microcentrifuge tube and add **100 µL PPS**. Mix by inverting the tube 20 times.

Note: Do not vortex the tube at this step

7. Centrifuge at 14,000 x g for 5 mins to collect pellet at the bottom of tube.

8. Transfer 400 µL supernatant carefully to well #2 or #8 of the pre-filled reagent plate.

9. Run the program (**Magbead_Feces**) on MPure-32[™] aNAP System. Alternatively, set a program on the instrument with the following settings:

| Step | Well | Process | Time (s) | | | Mixing Speed | Temp (°C) |
|------|--------|----------------------------|----------|------|---------|--------------|-----------|
| | | | Mix | Wait | Attract | | |
| 1 | #1/#7 | Magnetic Beads Preparation | 60 | 0 | 120 | Medium | RT |
| 2 | #2/#8 | Bind | 600 | 0 | 150 | Medium | RT |
| 3 | #3/#9 | Wash 1 | 180 | 0 | 120 | Medium | RT |
| 4 | #4/#10 | Wash 2 | 180 | 0 | 120 | Medium | RT |
| 5 | #5/#11 | Repeat Wash 2 | 180 | 0 | 150 | Medium | RT |
| 6 | #5/#11 | Dry | 0 | 600 | 0 | - | RT |
| 7 | #6/#12 | Elute | 300 | 0 | 150 | Medium | 55 |

10. Transfer eluted DNA into a clean 1.5 mL microcentrifuge tube. DNA is now ready for PCR and other downstream applications. Store the purified nucleic acid at -20°C for extended periods.

Note: If magnetic beads are observed in eluted DNA, centrifuge the eluent at 14,000 x g for 3-5 mins and transfer the clear supernatant to a clean 1.5 mL tube.

Lysis

Purify

Bind,
Wash,
Eluate

Ordering Information

Equipment

Automation extraction system from low to high throughput

| Catalog No. | Product Name | Throughput |
|-------------|-----------------------|------------------|
| 07EMC043 | MPure-32™ aNAP System | Up to 32 samples |
| 07EMC044 | MPure-96™ aNAP System | Up to 96 samples |



Instruments for lysing and homogenizing soil samples

| Catalog No. | Product Name |
|-------------|-------------------|
| 116004500 | FastPrep™ Classic |
| 116005500 | FastPrep-24™ 5G |
| 116010500 | FastPrep-96™ |
| 116012500 | SuperFastPrep-2™ |

Reagents

Wide range of reagent kits for extracting and purifying various types of biological samples for downstream applications.



| Catalog No. | Product Name | Pack Size |
|-------------|--|-----------|
| 117033100 | (MPure-32) MagBeads FastDNA® Kit for Soil (Ready-to-Use) | 96 preps |
| 117033200 | (MPure-32) MagBeads FastDNA® Kit for Feces (Ready-to-Use) | 96 preps |
| 117033300 | (MPure-32) MagBeads FastDNA/RNA Kit for Virus (Ready-to-Use) | 96 preps |
| 117033400 | (MPure-32) MagBeads FastRNA Kit (Ready-to-Use) | 96 preps |
| 117033500 | (MPure-32) MagBeads FastRNA Kit for FFPE (Ready-to-Use) | 96 preps |
| 117033600 | (MPure-32) MagBeads FastDNA® Kit (Ready-to-Use) | 96 preps |
| 117033700 | (MPure-32) MagBeads FastDNA Kit for Blood (Ready-to-Use) | 96 preps |
| 117033800 | (MPure-32) MagBeads FastDNA Kit for FFPE (Ready-to-Use) | 96 preps |
| 117033900 | (MPure-32) MagBeads Fast Circulating DNA Kit (Ready-to-Use) | 96 preps |



MP BIOMEDICALS

APAC: +65 6775 0008 | custserv.ap@mpbio.com

EUROPE: 00800 7777 9999 | custserv.eur@mpbio.com

AMERICAS: 800 854 0530 | custserv.na@mpbio.com

Learn more at www.mpbio.com

