

HEV ELISA 4.0

Instructions for Use



Date of Revision: 2016-02
MBE0011-ENG-3

Note Changes Highlighted

REF 23540-096: (96 tests)

NAME AND INTENDED USE

The **MP Diagnostics HEV ELISA 4.0** is an enzyme-linked immunosorbent assay intended for the detection of total antibodies to Hepatitis E Virus in human serum or plasma. It is intended as a screening test, requiring repeat testing of initially reactive specimens.

INTRODUCTION

Hepatitis E Virus (HEV) is a single-stranded, positive sense, non-enveloped RNA virus, which was first identified as an enterically-transmitted non-A, non-B hepatitis virus by Genelabs in 1990 (1,2). The course of the infection of HEV is generally acute and self-limiting without chronic sequelae. There is, however, a high incidence of mortality in pregnant women in the third trimester, about 10-20% (3) and a mortality rate of 1-2% in the general population, which is 10 times that of hepatitis A (HAV). With the cloning of the etiological agent of ET-NANBH at Genelabs and the identification of type common viral epitopes (1,2), specific diagnostic tools have been developed to detect antibodies to HEV.

Major epidemics of enterically transmitted non-A, non-B hepatitis (ET-NANBH) have been found to occur in developing regions such as Asia, the former USSR, Central America and Africa (3,4). Sporadic cases have been reported in developed nations, including Australia, the United Kingdom and the United States (5,6,7). Cases in developed nations have generally been associated with travel to endemic regions. However, accumulated evidences suggest that sporadic cases of HEV infections without an association with endemic regions also occur in a wide range of non-endemic areas, including Western Europe, Greece, United States, Australia, and Taiwan (8-17).

It has been demonstrated in the experiments that human HEV is capable of infecting animal species (18-21), while non-human primate may get infected with swine HEV (21). Recent studies on prevalence of HEV infection in animal show the high seroprevalence of antibody to HEV in different animal species, including swine, equine, roden, etc. Mounting evidences indicate that wide spread of HEV infection in animals, in particular swines, could represent an important reservoir for virus transmission. Some of the sporadic cases of HEV infection in non-endemic areas may be attributed to zoonotic transmission.

A HEV ELISA that is highly sensitive and specific is needed for the detection of HEV total antibodies in human serum or plasma.

The **MP Diagnostics HEV ELISA 4.0** utilises a proprietary recombinant antigen, which is highly conserved between different HEV strains (22,23,24), to detect the presence of specific antibodies including IgG, IgM and IgA against HEV.

DESCRIPTION OF SYMBOLS USED

The following are graphical symbols used in or found on MP Diagnostics products and packaging. These symbols are the most common ones appearing on medical devices and their packaging. Some of the common symbols are explained in more detail in the European and International Standard EN ISO 15223: 2012.

	Use by		In vitro diagnostic medical device
	Batch Code		Catalogue Number
	Temperature Limitation		Caution
	Manufacturer		Authorized Representative in the European Community
	Sufficient for <n> tests		Consult Instructions for Use
	Do not reuse		
	Contents		

CHEMICAL & BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The wells of the polystyrene microplate strips are coated with a proprietary recombinant antigen presenting a conformational epitope that is highly conserved between different HEV strains. The HRP conjugate is produced with the same recombinant antigen labeled with horseradish peroxidase. This conjugate is first diluted appropriately in diluent buffer prior to being dispensed into the antigen-coated wells of the microplates. Serum or plasma samples are then added to the antigen-coated wells containing the diluent buffer and the conjugate. After incubation, HEV specific antibodies (IgG, IgM and IgA), if present, will bind to both the antigens immobilised on the wells and the antigen of the conjugate in the diluent. Subsequently, the wells are thoroughly washed to remove the unbound materials. A substrate solution containing 3,3',5,5'-tetramethylbenzidine (TMB) is then added to each well. The presence of specific antibodies is indicated by the presence of blue colour solution after incubation. Reaction is terminated by addition of sulphuric acid. The colour intensity of the resulting yellow reaction product is measured at 450nm using microplate reader and its corresponding optical density or absorbance is proportional to the amount of antibodies present in the specimen.

KIT COMPONENTS

Component Description	Quantity Provided
HEV MICROPLATE Twelve 8-well strips per plate, sealed in an aluminum pouch with desiccant. Each microplate well contains adsorbed recombinant HEV protein. Store at 2°C to 8°C.	1 plate (96 tests)
NON-REACTIVE CONTROL Inactivated normal human serum, non-reactive for anti-HCV, anti-HEV, HBsAg and anti-HIV-1. Contains thimerosal and sodium azide as preservatives. Store at 2°C to 8°C.	1 vial (400µl)
REACTIVE CONTROL Inactivated human serum containing a high titer of IgG antibodies specific for HEV. Contains thimerosal and sodium azide as preservatives. Store at 2°C to 8°C.	1 vial (400µl)
SAM DILUENT (SAM = Sample Addition Monitor) Tris based saline solution containing heat-treated normal goat serum, bovine serum albumin and stabilizers. Contains BRONIDOX®L as preservative. Store at 2°C to 8°C.	1 bottle (100ml)
PLATE WASH CONCENTRATE (20x) Phosphate buffered saline with Tween-20. Contains chloroacetamide as preservative. Store at 2°C to 8°C.	1 bottle (120ml)
CONJUGATE HEV antigen labeled with horseradish peroxidase. Contains 0.02% thimerosal as preservative. Store at 2°C to 8°C.	1 vial (50µl)
SUBSTRATE BUFFER Buffer containing 3,3',5,5'-tetramethylbenzidine (TMB). Store in the dark at 2°C to 8°C.	1 bottle (12.5ml)
STOP SOLUTION 2M sulphuric acid solution. Store at 2°C to 8°C.	1 bottle (30ml)
PLATE COVERS Adhesive covers for microplate during incubation.	4 pieces
INSTRUCTIONS FOR USE	1 copy

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- For Professional use only.
- Please refer to the product labelling for information on potentially hazardous components.

HEALTH AND SAFETY INFORMATION



CAUTION: This kit contains materials of human origin. No test method can offer complete assurance that human blood products will not transmit infection.

HANDLE ASSAY SPECIMENS, REACTIVE AND NON-REACTIVE CONTROLS AS POTENTIALLY INFECTIOUS AGENTS. It is recommended that the components and test specimens be handled using good laboratory working practices. They should be disposed of in accordance with established safety procedures.

The **Reactive Control** and **Non-Reactive Control** contain 0.005% Thimerosal and 0.1% Sodium Azide. Sodium Azide can react with copper and lead used in some plumbing systems to form explosive salts. The quantities used in this kit are small, nevertheless when disposing of azide-containing materials they should be flushed away with relatively large quantities of water to prevent metal azide buildup in plumbing system.

Pursuant to EC regulation 1272/2008 (CLP), hazardous components are classified and labelled as follows:

Component:	Plate Wash Concentrate (20x)
Signal Word:	Warning
Pictogram:	
Hazard Statements:	H317 May cause an allergic skin reaction
Precautionary Statements:	P261 Avoid breathing dust/fume/gas/mist/vapours/spray. P272 Contaminated work clothing should not be allowed out of the workplace. P302+P352 IF ON SKIN: Wash with plenty of soap and water. P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
Supplemental Statements:	EUH210 Safety Data Sheet is available on request
Contains:	2% Chloroacetamide

Component:	Stop solution
Signal Word:	Danger
Pictogram:	
Hazard Statements:	H315 Causes skin irritation. H319 Causes serious eye irritation.

Precautionary Statements:	P264 Wash hands thoroughly after handling. P280 Wear protective gloves/protective clothing/eye protection/face protection. P312 Call a POISON CENTER or doctor/physician if you feel unwell. P362 Take off contaminated clothing and wash before reuse. P302+P352 IF ON SKIN: Wash with plenty of soap and water. P332+P313 If skin irritation occurs: Get medical advice/attention. P337+P313 If eye irritation persists: Get medical advice/attention. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
Supplemental Statements:	EUH210 Safety Data Sheet is available on request
Contains:	11.2% Sulphuric acid

- Avoid microbial contamination of reagents when opening and removing aliquots from the original vials or bottles.
- Do not pipette with mouth.
- Handle assay specimens, microplates, Reactive and Non-Reactive Controls as potentially infectious agents.
- Wear laboratory coats and disposable gloves while performing the assay. Discard gloves in bio-hazard waste-bags. Wash hands thoroughly afterwards.
- It is highly recommended that this assay be performed in a biohazard cabinet.
- Keep materials away from food and drink.
- In case of an accident or contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Consult a physician immediately in the event that contaminated materials are ingested or come in contact with open lacerations, or other breaks in the skin.
- Sulphuric acid can cause burns. **AVOID CONTACT.** If it comes into contact with skin, wash thoroughly with water.
- Avoid contact of sulphuric acid with any oxidizing agent or metal.
- Do not expose substrate to strong light.
- Wipe spills of potentially infectious materials immediately with absorbent paper and swab the contaminated area with an effective disinfecting agent before work is resumed.

ANALYTICAL PRECAUTIONS

- Use only sera or plasma samples collected in EDTA, Heparin, Sodium Citrate, K-Oxalate or Acid Citrate Dextrose (ACD). Before storage, ensure that blood clot or blood cells have been separated by centrifugation.
- Do not use whole blood or other body fluids.

STORAGE

- Store **MP Diagnostics HEV ELISA 4.0** kit and its components at 2°C to 8°C when not in use.
- All test reagents and strips in the closed or unopened condition, when stored at 2°C to 8°C, are stable until the expiry date given on the kit. Do not freeze the reagents.
- Crystals may form when Plate Wash Concentrate (20x) is stored at 2°C to 8°C. These must be dissolved by warming at 37°C prior to use.
- Precipitate may form when the Diluent is stored at 2°C to 8°C. This will not affect the performance of the kit.

SPECIMEN COLLECTION, TRANSPORT AND STORAGE

Serum or plasma specimens collected in EDTA, heparin, sodium citrate, K-oxalate or ACD may be used. Before storage, ensure that blood clot or blood cells have been separated by centrifugation.

Fresh specimens are preferred, specimens that undergo freeze-thaw cycles repeatedly are not recommended. Specimens should be stored at 2°C to 8°C if the test is to be run within 7 days of collection or frozen at ≤ -20°C if the test is to be delayed for more than 7 days. In addition, up to 0.1% Sodium Azide may be used to stabilize serum or plasma specimens stored at 2°C to 8°C.

Clear, non-haemolysed samples are preferred. Lipemic, icteric or contaminated (particulate) samples should be filtered (0.45µm) or centrifuged before testing.

Samples can be inactivated but this is not a requirement for optimal test performance.

Inactivate as follows:

- Loosen cap of sample container.
- Heat-inactivate sample at 56°C for 30 minutes in a water bath.
- Allow sample to cool down before retightening cap.
- Sample can be stored frozen until analysis.

Repeated freeze-thawing of sample is not recommended.

ADDITIONAL MATERIALS REQUIRED BUT NOT PROVIDED

- Disposable absorbent bench top paper and paper towels.
- Polypropylene tubes or containers.
- Graduated pipettes: 5ml, 10ml.
- Multichannel pipettor capable of delivering 20µl, 100µl, and 200µl.
- Pipettor capable of delivering 1-1000µl.
- Disposable pipette tips.
- Reagent reservoirs (troughs) with a capacity of 25ml.
- Deionised or distilled water, reagent grade quality.
- Flasks: 500ml, 1 litre.
- ELISA Microplate Washer. Alternatively, washing can be performed manually by using a multichannel pipettor delivering 0.3ml volumes and an aspirator device.
- A 37 ± 1°C incubator.
- A dual (A₄₅₀-A₆₂₀) or single (A₄₅₀) wavelength microplate reader.
- Effective disinfecting agent.

PREPARATION OF REAGENTS

- WORKING CONJUGATE**
 - WORKING CONJUGATE** should be prepared fresh prior to use.
 - Mix **CONJUGATE** and **DILUENT** thoroughly before use. **DO NOT SPIN** the mixture.
 - Dilute **CONJUGATE** at 1:500 dilution factor with **DILUENT**. For example, add 6.0µl conjugate into 3.0ml diluent.
 - Use only polypropylene containers or tubes.
 - 9.0ml of **WORKING CONJUGATE** is required for one microplate.

CONJUGATE PREPARATION CHART (1:500 dilution factor)		
Number of tests	Vol. of Conjugate (µl)	Vol. of Diluent (ml)
24	6.0	3.0
48	10.0	5.0
72	14.0	7.0
96	18.0	9.0

- DILUTED WASH BUFFER**
 - DILUTED WASH BUFFER** should be prepared fresh prior to use.
 - Dilute 1 volume of **PLATE WASH CONCENTRATE** with 19 volumes of distilled water (reagent grade quality). Mix well. Approximately 200ml of wash buffer is required to wash 1 plate.

ASSAY PROCEDURE

IMPORTANT: - Immunoassays of this nature are temperature-sensitive and time-dependent. Strict adherence to the assay procedure will ensure optimal assay performance. Deviations from the recommended procedure may lead to aberrant results.

- Prepare **WORKING CONJUGATE** as described in the **PREPARATION OF REAGENTS**.
- Remove microplate from the aluminum bag.
- Shake specimen and control vials before use.
- Fill a reagent reservoir with **WORKING CONJUGATE**. Using a multichannel pipettor, add 80µl of **WORKING CONJUGATE** to all wells. 80µl
- Wells A1 and B1 are 'BLANKS'. **DO NOT ADD SPECIMEN TO THESE WELLS.** Add 20µl of diluent per well to these wells. 20µl
- Add 20µl of specimen to the assigned well, starting at well A2. This will give a final specimen dilution of 1:5. Mix by pipetting up and down once. **DO NOT PLACE SPECIMEN IN AN EMPTY WELL.** 20µl

