



MP Biomedicals
Diagnostics Division
29525 Fountain Parkway
Solon, OH 44139

Prolactin ChLIA Kit
Catalog No.: 07M775A, 07M775B

INTENDED USE

The MP Biomedicals Prolactin ChLIA Kit is intended to be used for the quantitative determination of Prolactin Hormone Concentration in Human Serum by a Microplate Chemiluminescence Assay (ChLIA). This test is for in vitro diagnostic use only.

SUMMARY AND EXPLANATION OF THE TEST

Prolactin hormone (PRL), secreted from the lactotrophs of the anterior pituitary, is a protein consisting of a single polypeptide chain containing approximately 200 amino acids. The primary biological action of the hormone is on the mammary gland where it is involved in the growth of the gland and in the induction and maintenance of milk production. There is evidence to suggest that prolactin may be involved in steroidogenesis in the gonad, acting synergistically with luteinizing hormone (LH). High levels of prolactin appear to inhibit steroidogenesis as well as inhibiting LH and follicle stimulating hormone (FSH) synthesis at the pituitary gland (1, 2).

The clinical usefulness of the measurement of prolactin hormone (PRL) in ascertaining the diagnosis of hyperprolactinemia and for the subsequent monitoring the effectiveness of the treatment has been well established (3, 4).

In this method, PRL calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated monoclonal and enzyme labeled antibodies (directed against distinct and different epitopes of PRL) are added and the reactants mixed. Reaction between the various PRL antibodies and native PRL forms a sandwich complex that binds with the streptavidin coated to the well.

After the completion of the required incubation period, the enzyme-prolactin hormone antibody bound conjugate is separated from the unbound enzyme-prolactin hormone conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce light.


The employment of several serum references of known prolactin hormone levels permits the construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with prolactin hormone concentration.

PRINCIPLE

Immunoenzymometric assay (TYPE 3):

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme labeled and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti-PRL antibody.

Upon mixing monoclonal biotinylated antibody, the enzyme-labeled antibody and a serum containing the native antigen, reaction results between the native antigen and the antibodies, without competition or steric hindrance, to form a soluble sandwich complex. The interaction is illustrated by the following equation:



$\text{B}^{\text{m}}\text{Ab}(\text{m})$ = Biotinylated Monoclonal Antibody (Excess Quantity)

Ag_{PRL} = Native Antigen (Variable Quantity)

$\text{Enz}_{\text{Ab}(\text{p})}$ = Enzyme labeled Antibody (Excess Quantity)

$\text{Enz}_{\text{Ab}(\text{p})}\text{-Ag}_{\text{PRL}}\text{-B}^{\text{m}}\text{Ab}(\text{m})$ = Antigen-Antibodies Sandwich Complex

k_a = Rate Constant of Association

k_{-a} = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:



$\text{Streptavidin}_{\text{c.w.}}$ = Streptavidin immobilized on well

Immobilized complex = sandwich complex bound to the well

After a suitable time, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity, determined by reaction with a substrate that generates light, in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

REAGENTS

Materials Provided:

A. PRL Calibrators –1 mL/vial

Six (6) vials of references for PRL antigen in serum at levels of 0(A), 5(B), 10(C), 25(D), 50(E) and 100(F) ng/mL. Store at 2-8°C. A preservative has been added.

Note: The calibrators, human serum based, were calibrated using a reference preparation, which was assayed against the WHO 3rd, IS (84/500).

B. hPRL Tracer Reagent —13 mL/vial

One (1) vial containing enzyme labeled antibody, biotinylated monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

C. Light Reaction Wells – 96 wells

One 96-well white microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

D. Wash Solution Concentrate – 20 mL

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

E. Signal Reagent A – 7 mL/vial

One (1) bottle containing luminol in buffer. Store at 2-8°C.

F. Signal Reagent B – 7 mL/vial

One (1) bottle containing hydrogen peroxide (H₂O₂) in buffer. Store at 2-8°C.

G. Product Instructions

Note 1: Do not use reagents beyond the kit expiration date.

Note 2: Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

Note 3: Above reagents are for a single 96-well microplate

4.1 Required But Not Provided:

1. Pipette capable of delivering 50 µL volumes with a precision of better than 1.5%.
2. Dispenser(s) for repetitive deliveries of 0.100 mL and 0.350 mL volumes with a precision of better than 1.5%.
3. Microplate washers or a squeeze bottle (optional).
4. Microplate luminometer.
5. Test tube(s) for mixing substrates A&B.
6. Absorbent Paper for blotting the microplate wells.
7. Plastic wrap or microplate cover for incubation steps.
8. Vacuum aspirator (optional) for wash steps.
9. Timer.
10. Quality control materials

PRECAUTIONS

For In Vitro Diagnostic Use
Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

Safe Disposal of kit components must be according to local regulatory and statutory requirements.

SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood, serum in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.100ml of the specimen is required.

REAGENT PREPARATION

1. **Wash Buffer**
Dilute contents of Wash Concentrate to 1000 mL with distilled or deionized water in a suitable storage container. Store diluted buffer at 2-30°C for up to 60 days.
2. **Working Signal Reagent Solution - Store at 2 - 30°C.**
Determine the amount of reagent needed and prepare by mixing equal portions of Signal Reagent A and Signal Reagent B in a clean container. For example, add 1 mL of A and 1 mL of B per two (2) eight well strips (A slight excess of solution is made). **Discard the unused portion if not used within 36 hours after mixing.** If complete utilization of the reagents is anticipated, within the above time constraint, pour the contents of Signal Reagent B into Signal Reagent A and label accordingly.

Note: Do not use reagents that are contaminated or have bacteria growth.

TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum references and controls to room temperature (20-27°C).

****Test Procedure should be performed by a skilled individual or trained professional****

1. Format the microplate wells for each serum reference, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C**
2. Pipette 0.025 mL (25 µL) of the appropriate serum reference, control or specimen into the assigned well.
3. Add 0.100 mL (100 µL) of hPRL Tracer Reagent to all wells.
4. Swirl the plate gently for 20-30 seconds to mix and cover.
5. Incubate 45 minutes at room temperature.
6. Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
7. Add 350 µL of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat four (4) additional times for a total of five (5) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat four (4) additional times.**
8. Add 0.100 mL (100 µL) of working signal reagent to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells**
9. Incubate at room temperature in the dark for five (5) min.
10. Read the relative light units in each well, for minimum 0.5 – 1.0 seconds, using a microplate luminometer. **The results should be read within thirty (30) minutes of adding the substrate solution.**

QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and elevated ranges for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control records should be maintained and used to monitor batch to batch consistency.

CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of prolactin (PRL) concentration in unknown specimens.


1. Record the RLUs (Relative Light Unit) obtained from the printout of the microplate luminometer as outlined in Example 1.
2. Plot the RLUs for each duplicate serum reference versus the corresponding PRL concentration in ng/ml on linear graph paper.
3. Draw the best-fit curve through the plotted points.
4. To determine the concentration of hPRL for an unknown, locate the average RLU's of the unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average RLU's (33555) of the unknown intersects the calibration curve at (13.9 ng/mL) PRL concentration (See Figure 1).

Note: Computer data reduction software designed for chemiluminescence assays may also be used for the data reduction. **If such software is utilized, the validation of the software should be ascertained.**

EXAMPLE 1

Sample I.D.	Well Number	RLU (A)	Mean RLU (B)	Value (ng/mL)
Cal A	A1	113	139	0
	B1	164		
Cal B	C1	12028	12485	5
	D1	12941		
Cal C	E1	24693	24603	10
	F1	24513		
Cal D	G1	53221	53344	25
	H1	53468		
Cal E	A2	76850	77335	50
	B2	77820		
Cal F	C2	98568	100000	100
	D2	101432		
Ctrl 1	E2	14555	14690	5.9
	F2	14825		
Ctrl 2	G2	42680	42186	18.3
	H2	41692		
Patient	A3	32955	33555	13.9
	B3	34155		

Figure 1



* The data presented in Example 1 and Figure 1 is for illustration only and should not be used in lieu of a dose response curve prepared with each assay. In addition, the RLU's of the calibrators have been normalized to 100,000 RLU's for the F calibrator (greatest light output). This conversion minimizes differences caused by efficiency of the various instruments that can be used to measure light output.

LIMITATIONS OF PROCEDURE

Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of signal reagent initiates a kinetic reaction; therefore the signal reagent(s) should be added in the same sequence to eliminate any time-deviation during reaction.

- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Patient specimens with abnormally high prolactin levels can cause a hook effect, that is, paradoxical low absorbance results. If this is suspected, dilute the specimen 1/100 with '0' calibrator; reassay (multiply the result by 100). However, values as high as 3000 ng/mL have been found to absorb greater than the absorbance of the highest calibrator
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from MP Biomedicals' IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.

Interpretation

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.**
- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **MP Biomedicals shall have no liability.**
- If computer-controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.
- Patients receiving preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA) and may show either falsely elevated or depressed values when assayed.
- Pregnancy, lactation, and the administration of oral contraceptives can cause an increase in the level of prolactin.
- Drugs such as morphine, reserpine and the psychotropic drugs increase prolactin secretion (5, 6, and 7).
- Since prolactin hormone concentration is dependent upon diverse factors other than pituitary homeostasis, the determination alone is not sufficient to assess clinical status.

EXPECTED RANGES OF VALUES

A study of an apparent normal adult population was undertaken to determine expected values for the Prolactin ChLIA test. The expected values (95% confidence intervals) are presented in Table 1.

TABLE 1 Expected Values for the Prolactin ChLIA (in ng/mL)	
	Women
Adult (Number = 70)	1.2 -- 19.5
Postmenopausal (Number = 10)	1.5 -- 18.5
	Men
Adult (Number = 50)	1.8 -- 17.0

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

PERFORMANCE CHARACTERISTICS

Precision

The within and between assay precision of the Prolactin ChLIA test were determined by analyses on three different levels of control sera. The number, mean value, standard deviation (σ) and coefficient of variation for each of these control sera are presented in Table 2 and Table 3.

**TABLE 2
Within Assay Precision (Values in ng/mL)**

Sample	N	X	σ	C.V.
Level 1	20	5.4	0.23	4.3%
Level 2	20	18.4	0.67	3.6%
Level 3	20	40.8	2.78	6.8%

**TABLE 3
Between Assay Precision* (Values in ng/mL)**

Sample	N	X	σ	C.V.
Level 1	20	5.8	0.57	9.8%
Level 2	20	19.8	1.73	8.8%
Level 3	20	43.8	2.97	6.8%

*As measured in ten experiments in duplicate.

Sensitivity

The sensitivity (detection limit) was ascertained by determining the variability of the 0ng/ml serum calibrator and using the 2 σ (95% certainty) statistic to calculate the minimum dose. It was determined to be 0.11 ng/mL.

Accuracy

The Prolactin ChLIA test was compared with a reference Elisa method. Biological specimens from normal and abnormal populations were assayed. The total number of such specimens was 85. The least square regression equation and the correlation coefficient were computed for the Prolactin ChLIA test comparison with the reference method. The data obtained is displayed in Table 4.

**TABLE 4
Least Square Regression Analysis**

Method	Mean (x)	Analysis	Correlation Coefficient
This Method	18.5	$y = -1.63 + 1.01(x)$	0.978
Reference	19.4		

Only slight amounts of bias between this procedure and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

Specificity

The cross-reactivity of the Prolactin ChLIA test to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of prolactin hormone needed to produce the same light intensity.

Substance	Cross Reactivity	Concentration
Prolactin Hormone (PRL)	1.0000	--
Luteinizing Hormone (LH)	< 0.0001	1000 ng/mL
Follitropin (FSH)	< 0.0001	1000 ng/mL
Chorionic gonadotropin (CG)	< 0.0001	1000 ng/mL
Thyrotropin (TSH)	< 0.0001	1000 ng/mL
Growth Hormone (GH)	< 0.0001	1000 ng/mL

REFERENCES

- Maddox, P.R., Jones, D.L., Mansel, R.E., *Acta Endocrinol*, **125**, 621 (1991).
- Gonzales, E.R., *JAMA* **242**, 401 (1979).
- Tolis, G., *Hosp. Pract.* **15**, 85 (1980).
- Balagura, S., Frantz, A.G., Houseplan, E.M., *J Neurosurg* **51**, 42 (1979).
- Friessen, H., Hwang, P., *Ann Rev Med*, **24**, 251 (1973).
- Frantz, A. G., *N Eng J Med*, **298**, 201 (1978).

- Parke, D.N., *J. Med.* **301**, 873. (1979)
- Kao, P.C., Jiang, N.C. and Abboud, C.F., 'Radioimmunoassay of human homologous prolactin in serum with commercially available reagents' *Clin Chem.* **23**, 1563-1568. (1977)
- Haus, E., Lakatua, D.J., Halberg, F., Halberg, E., Cornelissen, G., Sackett, L.L. et al.: 'Chronobiological studies of plasma prolactin in women in Kyushu, Japan and Minnesota, USA.' *J.Clin. Endocrinol Metab.* **51**, 632-640. (1980)
- Christensen, J.M., Poulsen, O.M. and Anglov, T., 'Method, evaluation, quality control, and extreme quality assurance systems of analytical procedures', in: Sellaer HG, Seigel, H Eds. *Handbook on Metals in Clinical and Analytical Chemistry*. New York: Marcel Dekker 45-61 (1994).
- Phillips, G.B.: 'Relationship between serum levels of dehydroepiandrosterone sulfate, androstenedione and sex hormone in men and women.' *Eur.J.Endocrinol* **134**, 201-206. (1996)
- Toutiou Y., Carayon, A., Reinberg, A., Bogdan, A., Beck, H., 'Differences in seasonal rhythmicity of plasma prolactin in elderly human subjects: Detection in women but not in men', *J Endocrinol*, **96**, 65-71 (1983).
- Costonga, J., Janson PCV, Hermans, J., Van Wersch, JWJ and Bombacher P.J.: 'Short term and long term intra-individual variations and critical differences of clinical laboratory parameters', *J.Clin.Chem.& Clin.Biochem.* **985**, 23, 7-16.
- John R., McDowell, IFW, Scanton. MF and Ellis, AR.: 'Macroprolactin reactivities in prolactin assays: an issue for clinical laboratories and equipment manufacturers', *Clin Chem*, **46**, 884-885 (2000).
- Lindstedt, G., 'Endogenous antibodies against prolactin – a "new" cause of hyperprolactinemia', *Eur. J. Endo* **130**, 439-42 (1994).
- Cavaco, B., Prazeres, S., Santos, M.A., Sobrinho, L.G. and Leite, V., 'Hyperprolactinemia due to big big prolactin is differently detected by commercially available immunoassays', *J Endocrinol Invest*, **22**, 203-208 (1999).

Contact: **MP Biomedicals, LLC**
Diagnostics Division
29525 Fountain Parkway
Solon, OH 44139

Customer Service: (800) 854-0530
Fax: (440) 337-1180
www.mpbio.com

Please visit our website for other products and information.



MP Biomedicals Germany GmbH
Thüringer Straße 15
37269 Eschwege
GERMANY
Tel: 0049 5651 921 0
Fax: 0049 5651 921 181

QER NO.: Q19-014
04/19
DCO: 0642-R4