

MagBeads FastDNA™ Kit for Soil

Cat. No.:116561050 (50 preps) & 116561005 (5 preps)



Quick-Start Protocol

Revision 1.0 Mar 2021

Notes before starting

- Store Magnetic Beads at 2-8 °C upon arrival; do not freeze.
- Expect precipitation in Lysis Buffer S1; warming the solution to 55 °C will help dissolve the precipitates.
- Add 35 mL isopropanol to Binding Buffer MS and mark on the bottle.
- Add 50 mL 100% ethanol to Wash Buffer S and mark on the bottle.
- Vortex the sample in a Lysing Matrix E tube at maximum speed for 10 mins if a FastPrep® Instrument is unavailable. Secure samples on the vortex through an adapter to ensure homogenization.
- Centrifugation speed stated in the manual will be a guideline; use the maximum speed available if 14,000 x g is not feasible.

Manual Extraction

Lyse	<ol style="list-style-type: none">1 Add 100-500 mg soil sample to a Lysing Matrix E tube. Note: After adding the sample to the tube, ensure there is still 1/3-1/4 empty space remaining in the tube.2 Add 980 μL Lysis Buffer S1, 120 μL Lysis Buffer S2 and 10 μL RNase A Solution to the sample in the Lysing Matrix E tube and vortex 5-10 seconds to mix.3 Homogenize sample in a FastPrep® Instrument for 20-40 seconds at a speed setting of 6.0 m/s. Note: The speed and time can be changed according to different soil samples. Vortex 5-10 mins at maximum speed if a FastPrep® Instrument is not available. If homogenizers from other manufacturers are used, consult the Instruction Manual or manufacturer for appropriate homogenization parameters.4 Centrifuge at 14,000 x g for 5 mins to pellet debris.
Purify	<ol style="list-style-type: none">5 Transfer the supernatant (~800 μL) to a clean 1.5 mL microcentrifuge tube. Add 250 μL Inhibitor Removal MS and mix by inverting the tube 20 times.6 Centrifuge at 14,000 x g for 5 mins to pellet precipitate.
Bind	<ol style="list-style-type: none">7 Transfer the supernatant (~800 μL) to a new 2.0 mL microcentrifuge tube. Add an equal volume of Binding Buffer MS and 5 μL Magnetic Beads to the supernatant. Vortex or invert the tube to mix. Note: Ensure Magnetic Beads are thoroughly mixed before transferring to the supernatant.8 Place the tube on a shaker for 5 mins to allow binding.9 Place the tube on a magnetic rack for 3-5 mins, allow Magnetic Beads to settle, then discard supernatant. Note: If the supernatant is too turbid or Magnetic Beads are attracting slowly, extend attraction time.
Wash	<ol style="list-style-type: none">10 Add 800 μL Wash Buffer S to the tube and place on the shaker for 3 mins.11 Place the tube on the magnetic rack for 1 min, allow Magnetic Beads to settle, then discard supernatant.12 Repeat step 10 to step 11 for a second wash step.13 Air dry Magnetic Beads for 5-10 mins at 55 °C by placing the tube on a heat block. Note: This is for removal of residual ethanol and ensure Magnetic Beads are completely dry.
Elute	<ol style="list-style-type: none">14 Add 100 μL DES Buffer to resuspend Magnetic Beads and incubate on a heat block at 55 °C for 5 mins.15 Place the tube on a magnetic rack for 3-5 mins until Magnetic Beads have settled, and transfer the supernatant (eluted DNA) to a clean 1.5 mL microcentrifuge tube. DNA is now ready for PCR and other downstream applications. Store at -20 °C for extended periods. Note: If the supernatant is too turbid or there are still Magnetic Beads remaining, please centrifuge at 14,000 x g for 3-5 mins and transfer the supernatant again.

MagBeads FastDNA™ Kit for Soil

Cat. No.:116561050 (50 preps) & 116561005 (5 preps)



Quick-Start Protocol

Revision 1.0 Mar 2021

Automated Extraction

Lyse	1 Add 100-500 mg soil sample to a Lysing Matrix E tube. Note: After adding the sample to the tube, ensure there is still 1/3- 1/4 empty space remaining in the tube.																																																																		
	2 Add 980 µL Lysis Buffer S1 , 120 µL Lysis Buffer S2 and 10 µL RNase A Solution to the sample in the Lysing Matrix E tube and vortex 5- 10 seconds to mix.																																																																		
Purify	3 Homogenize sample in a FastPrep® Instrument for 20- 40 seconds at a speed setting of 6.0 m/s. Note: The speed and time can be changed according to different soil samples. Vortex 5- 10 mins at maximum speed if a FastPrep® Instrument is not available. If homogenizers from other manufacturers are used, consult the Instruction Manual or manufacturer for appropriate homogenization parameters.																																																																		
	4 Centrifuge at 14,000 x g for 5 mins to pellet debris.																																																																		
Auto Bind, Wash & Elute	5 Transfer the supernatant (~800 µ L) to a clean 1.5 mL microcentrifuge tube. Add 250 µ L Inhibitor Removal MS and mix by inverting the tube 20 times.																																																																		
	6 Centrifuge at 14,000 x g for 5 mins to pellet precipitate.																																																																		
Auto Bind, Wash & Elute	7 Transfer 400 µ L supernatant to well 2 and well 3 of a 96- well plate. Add additional reagents into the respective wells as shown below.																																																																		
	<table border="1"><thead><tr><th>Well</th><th>Reagents</th><th>Volume (µL)</th></tr></thead><tbody><tr><td rowspan="2">1</td><td>Deionized water</td><td>800</td></tr><tr><td>Magnetic Beads</td><td>5</td></tr><tr><td rowspan="2">2</td><td>Sample supernatant</td><td>400</td></tr><tr><td>Binding Buffer MS</td><td>400</td></tr><tr><td rowspan="2">3</td><td>Sample supernatant</td><td>400</td></tr><tr><td>Binding Buffer MS</td><td>400</td></tr><tr><td>4</td><td>Wash Buffer S</td><td>800</td></tr><tr><td>5</td><td>Wash Buffer S</td><td>800</td></tr><tr><td>6</td><td>DES Buffer</td><td>100</td></tr></tbody></table>	Well	Reagents	Volume (µL)	1	Deionized water	800	Magnetic Beads	5	2	Sample supernatant	400	Binding Buffer MS	400	3	Sample supernatant	400	Binding Buffer MS	400	4	Wash Buffer S	800	5	Wash Buffer S	800	6	DES Buffer	100																																							
Well	Reagents	Volume (µL)																																																																	
1	Deionized water	800																																																																	
	Magnetic Beads	5																																																																	
2	Sample supernatant	400																																																																	
	Binding Buffer MS	400																																																																	
3	Sample supernatant	400																																																																	
	Binding Buffer MS	400																																																																	
4	Wash Buffer S	800																																																																	
5	Wash Buffer S	800																																																																	
6	DES Buffer	100																																																																	
Auto Bind, Wash & Elute	8 Run the instrument according to following settings.																																																																		
	<table border="1"><thead><tr><th rowspan="2">Step</th><th rowspan="2">Well</th><th rowspan="2">Process</th><th colspan="3">Time (s)</th><th rowspan="2">Mixing Speed</th><th rowspan="2">Temp (°C)</th></tr><tr><th>Mix</th><th>Wait</th><th>Attract</th></tr></thead><tbody><tr><td>1</td><td>1</td><td>Magnetic Beads Preparation</td><td>60</td><td>0</td><td>120</td><td>Medium</td><td>RT</td></tr><tr><td>2</td><td>2</td><td>Bind</td><td>300</td><td>0</td><td>120</td><td>Medium</td><td>RT</td></tr><tr><td>3</td><td>3</td><td>Bind</td><td>300</td><td>0</td><td>120</td><td>Medium</td><td>RT</td></tr><tr><td>4</td><td>4</td><td>Wash 1</td><td>180</td><td>0</td><td>120</td><td>Medium</td><td>RT</td></tr><tr><td>5</td><td>5</td><td>Wash 2</td><td>180</td><td>0</td><td>120</td><td>Medium</td><td>RT</td></tr><tr><td>6</td><td>5</td><td>Dry</td><td>0</td><td>600</td><td>0</td><td>-</td><td>RT</td></tr><tr><td>7</td><td>6</td><td>Elute</td><td>300</td><td>0</td><td>150</td><td>Medium</td><td>55</td></tr></tbody></table>	Step	Well	Process	Time (s)			Mixing Speed	Temp (°C)	Mix	Wait	Attract	1	1	Magnetic Beads Preparation	60	0	120	Medium	RT	2	2	Bind	300	0	120	Medium	RT	3	3	Bind	300	0	120	Medium	RT	4	4	Wash 1	180	0	120	Medium	RT	5	5	Wash 2	180	0	120	Medium	RT	6	5	Dry	0	600	0	-	RT	7	6	Elute	300	0	150	Medium
Step	Well				Process	Time (s)				Mixing Speed	Temp (°C)																																																								
		Mix	Wait	Attract																																																															
1	1	Magnetic Beads Preparation	60	0	120	Medium	RT																																																												
2	2	Bind	300	0	120	Medium	RT																																																												
3	3	Bind	300	0	120	Medium	RT																																																												
4	4	Wash 1	180	0	120	Medium	RT																																																												
5	5	Wash 2	180	0	120	Medium	RT																																																												
6	5	Dry	0	600	0	-	RT																																																												
7	6	Elute	300	0	150	Medium	55																																																												
Auto Bind, Wash & Elute	9 Transfer eluted DNA from well 6 into a clean 1.5 mL microcentrifuge tube. DNA is now ready for PCR and other downstream applications. Store at -20 °C for extended periods.																																																																		
	Note: If there are still Magnetic Beads remaining in eluted DNA, please centrifuge at 14,000 x g for 3- 5 mins and transfer the supernatant again.																																																																		

Order Information

Product	Package	Cat. No.
MagBeads FastDNA™ Kit for Soil	50preps	116561050
FastPrep - 24™ 5G Instrument Including 24 x 2 mL samples adapter	1ea	116005500
MP Magnetic Rack 24	1ea	116570413
MP Magnetic Rack 8	1ea	116570426

