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	Private Label Manufacturer: MP Biomedicals Germany GmbH Thueringer Str. 15 37269 Eschwege, Germany	Original Equipment Manufacturer: Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimel North Industrial Park, Xiamen, Fujian, 361021, P.R. China
	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzhong Street, Jinzhou District, Dalian, 116100 Liaoning, China	
	Jiangsu Hanheng Medical Technology Co.,Ltd. 16-B4,#1 North Qingyang Road, Tianning District, 213017, Changzhou, Jiangsu, China	

MP Biomedicals Germany GmbH

laientest@mpbio.com

0049.5651.921111

Learn more at www.mpbio.com

English

RAPID SARS-COV-2 ANTIGEN TEST CARD

SELF-TEST Instructions for use

LAY TEST FOR ANTERIOR NASAL SWAB SPECIMENS

Catalog Number **REF**: 07AG6001BS (1 Test/Box); 07AG6005BS (5 Tests/Box); 07AG6020BS (20 Tests/Box)

- For private use/home use/self-testing.
- Please follow the instructions for use carefully.
- The repacking of 5 and 20 packs into smaller units is not permitted (separating is prohibited).

INTENDED USE

Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 viral antigen in anterior nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset. Rapid SARS-CoV-2 Antigen Test Card shall not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be assisted by an adult.

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection. Asymptomatic infected people can also be an infectious source. Based on current epidemiological studies, the incubation period is 1 to 14 days, but mostly 3 to 7 days. The main manifestations include fever, fatigue, loss of smell and/or taste, and a dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

MATERIALS PROVIDED

Components	For 1 Test/Box	For 5 Tests/Box	For 20 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5	20
Sterile swab	1	5	20
Extraction tube	1	5	20
Sample extraction buffer	1	5	20
Instructions for use (this leaflet)	1	1	1
Tube stand	Back of box	1	1

IMPORTANT INFORMATION BEFORE THE EXECUTION

- Read this instruction guide carefully.
- Do not use the product beyond the expiration date.
- Do not use the product if the pouch is damaged or the seal is broken.
- Store the test device at 4 to 30°C in the original sealed pouch. Do Not Freeze.
- The product should be used at room temperature (15°C to 30°C). If the product has been stored in a cool area (less than 15°C), leave it at normal room temperature for 30 minutes before using.
- Handle all specimens as potentially infectious.
- Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
- Use the swabs included in the test kit to ensure optimal performance of the test.
- Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.
- Blow the nose several times before collecting specimen.
- The specimens should be tested as soon as possible after collection.
- Apply the drops of test specimen only to the specimen well (S).
- Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
- Children under 14 years of age should be assisted by an adult.

LIMITATIONS

- The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as part of this test.
- Proper specimen collection is critical. Failure to follow the procedure may result in inaccurate test results.
- Improper collection, storage, or even freezing and thawing of the specimen can lead to inaccurate test results.
- If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after evaluation of all clinical and laboratory results.
- A negative result does not exclude viral infection except for SARS-CoV-2 and should be confirmed by molecular diagnostic methods (e.g. PCR) if COVID-19 is suspected.
- A positive result does not rule out an additional infection with other disease-causing agents.
- The SARS-CoV-2 rapid antigen test can detect both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and may not correlate with other diagnostic methods performed on the same specimen.
- Users should test specimens as soon as possible after specimen collection and within two hours of specimen collection.
- Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. The method of nasopharyngeal swab sampling should only be performed by healthcare professionals.
- It could be possible that virus mutations might be detected with lower sensitivity or not at all.
- The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.
- The kit was validated with the swabs provided. Use of alternative swabs may result in false negative results.
- The validity of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
- Cross-reactivity of the test was evaluated by testing viruses and other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Cross-Reactivity-Study. The listed viruses and other microorganisms have no effect on the test results, except the Human SARS-coronavirus. Positive test results do not rule out co-infections with other disease-causing agents. Positive results may occur in cases of infection with SARS-CoV.

PREPARATION

- Wash your hands.
- Clear clean and dry a flat surface.
- Check the kit contents.
- Make sure that nothing is damaged or broken.
- Have a timer ready.
- Blow your nose several times before taking the sample.
- Wash your hands again.

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Peel open the wrapper of the swab and carefully pull out the swab.

CAUTION: Do not touch the soft, fabric tip of the swab with your hands.

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Carefully insert the swab into one nostril. The swab tip should be inserted no less than 2.5 cm from the edge of the nostril. Roll the swab 3-4 times along the mucosa inside the nostril to ensure that both mucus and cells are collected. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril.

CAUTION: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.

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Place the swab with the sample into the extraction tube. Rotate the swab 3-5 times. **Leave the swab in the extraction solution for 1 minute.**

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Pinch the extraction tube together with fingers while removing the swab to leave as much

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	Private Label Manufacturer: MP Biomedicals Germany GmbH Thueninger Str. 15 37269 Eschwege, Allemagne		Original Equipment Manufacturer: Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian, 361021, P.R. China
	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yongzheng Street, Jinzhou District, Dalian, 116100 Liaoning, Chine		
	Jiangsu Hanheng Medical Technology Co., Ltd. 16-B4, #1 North Qingyang Road, Tanning District, 213017, Changzhou, Jiangsu, Chine		
			Luxus Lebenswelt GmbH Kochstr.1, 47877, Willich, Allemagne

MP Biomedicals Germany GmbH
laientest@mpbio.com
 0049.5651.921111

Learn more at www.mpbio.com