

# SPINeasy® DNA Kit for Yeast

For the isolation of DNA from Yeast sample

**Size: 50 and 5 PREPS**

**Storage: 15-25 °C**

**Cat. No.: 116557050 (50 PREPS) / 116557000 (5 PREPS)**

**Content Version: August 2023**

# Table of Contents

1. Introduction to SPINeasy® DNA Kit for Yeast.....	3
2. Kit Components and User Supplied Materials .....	4
3. Storage and Kit Stability .....	5
4. Important Consideration Before Use .....	5
5. Safety Precautions .....	6
6. Protocol.....	6
7. Flow Chart.....	8
8. Data .....	9
9. Troubleshooting .....	11
10. Product Use Limitation & Warranty .....	12

## 1. Introduction to SPINeasy® DNA Kit for Yeast

SPINeasy® DNA Kit for Yeast is a high-performance genomic DNA (gDNA) extraction kit which is developed based on silica-membrane spin-column technology. This kit enables quick isolation of gDNA from yeast cells, typically in less than 30 minutes. Provided in the kit, Column S and kit buffers are designed to deliver gDNA extracts of high yield and purity; the extracted gDNA is compatible with downstream applications such as PCR amplification, restriction enzyme digestion, and sequencing.

Visit [www.mpbio.com](http://www.mpbio.com) to explore additional products to support your research.

### Kit Specifications at a Glance

Technology	Silica membrane technology
Format	Mini spin column
Sample	Yeast
Sample amount	up to $4 \times 10^8$ cells.
Observed yield	up to 12 µg (number of cells and processing dependent)
Elution volume	50-200 µL
Preparation time	30 min

## 2. Kit Components and User Supplied Materials

### 2.1 SPINeasy® DNA Kit for Yeast Components

Product	50 reactions (Cat.No.: 116557050)		5 reactions (Cat.No.: 116557000)	
	Package	Cat. No.	Package	Cat. No.
Lysing Matrix YG	50 x 2 mL	116557051	5 x 2 mL	116557001
Lysis Buffer Y	30 mL	116557052	3.0 mL	116557002
Binding Buffer Y	30 mL	116557053	3.0 mL	116557003
Wash Buffer Y1	32 mL	116557054	3.2 mL	116557004
Wash Buffer Y2	4.8 mL	116557055	480 µL	116557005
Elution Buffer EB	12 mL	116552054	1.2 mL	116552004
RNase A	280 µL	116557056	28 µL	116557006
Column S	50 ea	116530058	5 ea	116530008
Collection tube	50 ea	116546059	5 ea	116558009
Elution tube	50 ea	116546060	5 ea	116546010
Quick-Start Protocol	1 each		1 each	
Instruction Manual	Available <a href="http://www.mpbio.com">www.mpbio.com</a>		Available <a href="http://www.mpbio.com">www.mpbio.com</a>	
MSDS & CoA	Available <a href="http://www.mpbio.com">www.mpbio.com</a>		Available <a href="http://www.mpbio.com">www.mpbio.com</a>	

### 2.2 User Supplied Materials

- FastPrep® Instrument - FastPrep-24™ 5G (Cat. No.116005500) or Vortex
- Microcentrifuge capable of centrifugation of at least 15,000 g
- Absolute ethanol
- Single-channel pipettors (5 µL-1000 µL)
- Nuclease-free, aerosol-preventive tips

### 3. Storage and Kit Stability

The SPINeasy® DNA Kit for Yeast is guaranteed until the expiry date stated on the kit when stored at room temperature (15-25°C). For extended storage or storage in dry condition (humidity < 40%), store the Column S at 2-8°C to maintain their performance.

### 4. Important Consideration Before Use

Please check as appropriate:

Add **40 mL** (50 PREPS kit) or **4 mL** (5 PREPS kit) of absolute ethanol to **Wash Buffer Y2**.

Mark on the bottle after addition.

Centrifugation speed stated in the manual will be a guideline; use the maximum speed available if 15,000 g is not feasible.

## 5. Safety Precautions

Lysis Buffer Y, Binding Buffer Y and Wash buffer Y1 contain a component that can be harmful if swallowed and may cause irritation when in contact with skin and eyes. To prevent accidental ingestion, do not eat, drink, or smoke when using this product. Wear personal protective equipment (gloves, lab coat and eye protection) to prevent contact with the skin or mucous membranes. Consult the Material Safety Data Sheet at [www.mpbio.com](http://www.mpbio.com) for additional details.

## 6. Protocol

### Things to do before starting.

- Add indicated amount of absolute ethanol to **Wash Buffer Y2**.

### Procedure:

#### 1. Sample Preparation

##### a. Harvesting cells from liquid cultures

Transfer the culture **1.0-2.0 mL (maximum  $4 \times 10^8$  cells)** to a microcentrifuge tube. Centrifuge for **5 min** at **15,000 g** and discard the supernatant.

##### b. Harvesting colonies from solid medium

Use an inoculating loop to transfer single colony ( $\geq 2$  mm in diameter) or colonies (amount  $< 50$  mg) from the solid medium.

2. Add **560  $\mu$ L Lysis Buffer Y** and **5  $\mu$ L RNase A** to each microcentrifuge tube of pellet. Resuspend pellet with the lysis buffer and RNase A mixture by pipetting repeatedly. Transfer all lysate to **Lysing matrix YG** tube.

3. Homogenize either in a **FastPrep Instrument** for **35 sec** at speed setting of **5.0 m/s** or vortex for **10 min** at speed setting of **2500 rpm**.

**Note:** For cells  $> 4 \times 10^8$  or colonies  $> 50$  mg, homogenize twice.

4. Centrifuge at **15,000 g** for **2 min**.
5. Slowly add **560  $\mu$ L Binding Buffer Y** to the mixture. Mix thoroughly by inversion.
6. Centrifuge at **15,000 g** for **1 min**.


**Note:** To achieve the best performance, add **200  $\mu$ L 1 M NaOH** (prepared by users) to the

*Column S before step 7, centrifuge at 15,000 g for 1 min. Discard the flow-through.*

7. Carefully transfer **750 µL** of the supernatant to **Column S** placed in a **Collection Tube** (provided). Centrifuge at **15,000 g** for **1 min**. Discard the flow-through.
8. Add **600 µL Wash Buffer Y1** to the wall of **Column S**. Centrifuge at **15,000 g** for **1 min**. Discard the flow-through.
9. Add **700 µL Wash Buffer Y2** to the wall of **Column S**. Centrifuge at **15,000 g** for **1 min**. Discard the flow-through.
10. Centrifuge at **15,000 g** for another **2 min** to dry the column.
11. Transfer **Column S** into a clean **1.5 mL Elution tube** (provided). Add **50-200 µL Elution Buffer EB** or ultrapure water (**pH > 7.0**) (not provided) to the center of **Column S**, incubate for **2 min**, and centrifuge at **15,000 g** for **1 min** to collect the eluted DNA.

**Optional:** *For maximum DNA recovery, reload the eluted DNA back to the center of **Column S** and centrifuge at **15,000 g** for **1 min** to collect the eluted DNA.*

## 7. Flow Chart






## 8. Data

SPINeasy® DNA Kit for Yeast delivers DNA yield with excellent purity that is in proportion with the input sample amount of *Saccharomyces cerevisiae* ( $1 \times 10^7$  to  $4 \times 10^8$  cells). The 260/280 and 260/230 ratios obtained from the extracted DNA are within the acceptable range (Figure 1).


SPINeasy® DNA Kit for Yeast delivers DNA yield doubled that obtained from Competitor Z. Despite Competitor Q and T giving higher spectrophotometric DNA yield, agarose gel electrophoresis result shows that the DNA samples are degraded. SPINeasy® DNA for Yeast gave higher purity ratios of 260/280 and 260/230 than all the competitor's kits (Figure 2).

**Figure 1:** SPINeasy® DNA Kit for Yeast generates DNA yield that is proportionate to the input sample amount, *Saccharomyces cerevisiae* ( $1 \times 10^7$  to  $4 \times 10^8$  cells). 3  $\mu$ L of purified genomic DNA from a 100  $\mu$ L eluate were analyzed by electrophoresis on a 1% agarose gel.



Sample	No.	Concentration (ng/ $\mu$ L)	A260/A280	A260/A230
1-2 x $10^7$ cells	1	12.784	1.87	2.16
	2	15.396	1.86	2.08
2-5 x $10^7$ cells	3	29.149	1.89	2.20
	4	27.321	1.87	2.17
1-2 x $10^8$ cells	5	66.683	1.88	2.27
	6	66.645	1.89	2.25
2-4 x $10^8$ cells	7	115.567	1.88	2.15
	8	125.671	1.90	2.06

**Figure 2: Genomic DNA extracted from *Saccharomyces cerevisiae* ( $5 \times 10^7$  cells).** The SPINeasy® DNA Kit for Yeast was used to extract genomic DNA (gDNA) from *Saccharomyces cerevisiae*. The yield and purity were assessed and compared with three different competitor kits (Q, T, Z). The MP Kit exhibited a purity range exceeding 1.8 while maintaining a high yield. Competitors Q and T yielded high quantities of gDNA but with a lower A260/A230 ratio. Competitor Z yielded both a lower quantity of gDNA and a lower purity ratio. Gel electrophoresis of the extracted gDNA displays maximum band integrity for the MP Kit.



Sample	No.	Concentration (ng/μL)	A260/A280	A260/A230
MP	1	60.363	1.87	1.91
	2	61.138	1.87	2.05
Competitor Q	3	72.729	2.05	0.62
	4	74.359	2.10	0.61
Competitor T	5	100.989	2.12	1.39
	6	95.493	2.14	1.61
Competitor Z	7	32.146	1.85	1.34
	8	22.194	1.75	1.02

## 9. Troubleshooting

Problem	Possible Cause	Recommendation
Low DNA Yield	Incomplete sample lysis or homogenization	Reduce the number of samples or increase the lysis time.
	DNA inefficiently eluted	Perform elution twice, 100 µL elution buffer each time, for a total elution volume of 200 µL.
	Water of incorrect pH used for elution	The low pH of deionized water from some water purifiers may reduce DNA yield. If eluting with water, ensure that the pH of the water is at least 7.0.
Smear DNA band	DNA not efficiently bound to column	Before use, equilibrate the column with 200 µL of 1 M NaOH, and centrifuge 1 minute at the maximum speed. Discard the flow-through.
	Yeast sample is too old	Use freshly collected yeast culture. Avoid repeated freeze/thaw cycles of the sample.
Low A260/A230 ratios	Mechanical sample disruption is too vigorous	If extracting from yeast cells $< 1 \times 10^7$ , alternative lysis method can be used: lysis by vortexing at the speed of 2000 rpm for 5 minutes.
	Carry over ethanol in eluant	Lengthen the column-drying time or increase the centrifugation speed to ensure complete drying of the column membrane and transfer the column into the elution tube carefully so that the column does not come into contact with the flow-through.
High A260/A230 ratios	Incomplete purification (inadequate lysis / overloading of column)	Reduce the amount of starting material
	EDTA in eluant	EDTA present in Elution Buffer EB may be associated with higher than usual A260/A230 ratio and this will not have any negative effect on downstream applications.
Poor PCR Performance	Too much DNA used	Check and adjust the DNA concentration.
	High inhibitor	Dilute template DNA to reduce the inhibitors.
	PCR protocol or reagents were not optimal	Verify PCR reagents and protocol with positive control; adjustment on reaction conditions or primer selection may be necessary.

## 10. Product Use Limitation & Warranty

The products presented in this instruction manual are for research or manufacturing use only. They are not to be used as drugs or medical devices to diagnose, cure, mitigate, treat or prevent diseases in humans or animals, either as part of an accepted course of therapy or in experimental clinical investigation. These products are not to be used as food, food additives or general household items. Purchase of MP Biomedicals products does not grant rights to reproduce, modify, or repackage the products or any derivative thereof to third parties. MP Biomedicals makes no warranty of any kind, expressed or implied, including merchantability or fitness for any particular purpose, except that the products sold will meet our specifications at the time of delivery.

Buyer's exclusive remedy and the sole liability of MP Biomedicals hereunder shall be limited to, at our discretion, no replacement or compensation, product credits, refund of the purchase price of, or the replacement of materials that do not meet our specification. By acceptance of the product, Buyer indemnifies and holds MP Biomedicals harmless against, and assumes all liability for, the consequence of its use or misuse by the Buyer, its employees, or others, including, but not limited to, the cost of handling. Said refund or replacement is conditioned on Buyer notifying within thirty (30) days of receipt of product. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

**Australia**

Tel: +61 2.8824.2100  
 Tel: +61 1800.249.998  
 Email: [custserv.au@mpbio.com](mailto:custserv.au@mpbio.com)

**Austria & Germany**

Tel: 0800.426.67.337  
 Tel: 00800.7777.9999  
 Email: [custserv.de@mpbio.com](mailto:custserv.de@mpbio.com)

**Belgium**

Tel: 00800.7777.9999  
 Email: [custserv.be@mpbio.com](mailto:custserv.be@mpbio.com)

**Canada**

Tel: +1 800.854.0530  
 Email: [custserv.ca@mpbio.com](mailto:custserv.ca@mpbio.com)

**China**

Tel: +86 400.150.0680  
 Email: [custserv.cn@mpbio.com](mailto:custserv.cn@mpbio.com)

**Europe**

Tel: +33 3.88.67.54.25  
 Tel: +33 00800.7777.9999  
 Email: [custserv.eur@mpbio.com](mailto:custserv.eur@mpbio.com)

**France**

Tel: +33 3.88.67.54.25  
 Email: [custserv.fr@mpbio.com](mailto:custserv.fr@mpbio.com)

**India**

Tel: +91 22.27636921/22/25  
 Email: [custserv.in@mpbio.com](mailto:custserv.in@mpbio.com)

**Italy**

Tel: 00800.7777.9999  
 Email: [custserv.it@mpbio.com](mailto:custserv.it@mpbio.com)

**Japan**

Tel: +81 3.6667.0730  
 Email: [custserv.jp@mpbio.com](mailto:custserv.jp@mpbio.com)

**Latin America**

Tel: +1 800.854.0530  
 Tel: +1 440.337.1200  
 Email: [custserv.la@mpbio.com](mailto:custserv.la@mpbio.com)

**New Zealand**

Tel: +64 9.912.2460  
 Email: [custserv.nz@mpbio.com](mailto:custserv.nz@mpbio.com)

**North America**

Tel: +1 800.854.0530  
 Tel: +1 440.337.1200  
 Email: [custserv.na@mpbio.com](mailto:custserv.na@mpbio.com)

**Poland**

Tel: 00800.7777.9999  
 Email: [custserv.po@mpbio.com](mailto:custserv.po@mpbio.com)

**Russia**

Tel: +7 495 604.13.44  
 Email: [custserv.rs@mpbio.com](mailto:custserv.rs@mpbio.com)

**Serbia**

Tel: +381 11.242.1972  
 Email: [custserv.se@mpbio.com](mailto:custserv.se@mpbio.com)

**Singapore/ APAC**

Tel: +65 6775.0008  
 Tel: +65 6394.7675  
 Email: [custserv.ap@mpbio.com](mailto:custserv.ap@mpbio.com)

**South Korea**

Tel: +82 2.425.5991  
 Email: [custserv.kr@mpbio.com](mailto:custserv.kr@mpbio.com)

**Switzerland**

Tel: 00800.7777.9999  
 Email: [custserv.ch@mpbio.com](mailto:custserv.ch@mpbio.com)

**The Netherlands**

Tel: 00800.7777.9999  
 Email: [custserv.nl@mpbio.com](mailto:custserv.nl@mpbio.com)

**United Kingdom**

Tel: 0800.282.474  
 Email: [custserv.uk@mpbio.com](mailto:custserv.uk@mpbio.com)



[www.mpbio.com](http://www.mpbio.com)