

SPINeasy RNA Kit for Tissue (With Lysing Matrix)

Spin Column Purification for Easy Isolation of RNA from
Animal Tissue, Plant Tissue, and Cell Culture Samples

Size: 50 & 5 preps

Storage: 15-25 °C

Cat. No.: 116543050 (50 PREPS)
116543000 (5 PREPS)

Content Version: November 2023



Table of Contents

| | |
|---|----|
| 1. Introduction to SPINeasy RNA Kit for Tissue (With Lysing Matrix) | 3 |
| 2. Kit Components and User Supplied Materials | 4 |
| 3. Storage and Kit Stability | 5 |
| 4. Important Consideration Before Use | 5 |
| 5. Safety Precautions | 6 |
| 6. Protocol..... | 7 |
| 7. Data | 9 |
| 8. Troubleshooting | 11 |
| 9. Product Use Limitation & Warranty | 12 |

1. Introduction to SPINeasy RNA Kit for Tissue (With Lysing Matrix)

SPINeasy RNA Kit for Tissue (With Lysing Matrix) is a silica-membrane spin-column kit that enables quick and convenient purification of total RNA from various tissue types from animals, plants and tissue cultures, without the use of toxic substances such as phenol and chloroform. The use of our specially formulated Lysis Buffer TR and Lysing Matrix A in combination with FastPrep® Instruments from MP Biomedicals enables highly efficient lysis of tissue samples within seconds. With a simple workflow, this kit allows multiple samples to be processed simultaneously. Total RNA of high quality and integrity can be purified typically within 40 minutes and is immediately available for downstream applications such as RT-PCR.

Visit www.mpbio.com to explore additional products to support your research.

Kit Specifications at a Glance

| | |
|------------------|--|
| Technology | Silica membrane technology |
| Format | Mini spin column |
| Sample | Animal tissue, plant tissue, and cell culture samples |
| Sample amount | 10-30 mg (animal tissue), 50-100 mg (plant tissue), 1x10 ⁶ cells (cell culture samples) |
| Elution volume | 100 µL |
| Preparation time | 40 min |

2. Kit Components and User Supplied Materials

2.1 SPINeasy RNA Kit for Tissue (With Lysing Matrix) Component

| Components | 50 PREPS (Cat.No.: 116543050) | | 5 PREPS (Cat.No.: 116543000) | |
|-------------------------------|--|-----------|---------------------------------|-----------|
| | Package | Cat. No. | Package | Cat. No. |
| Equilibration Buffer R | 12 mL | 112075704 | 1.2 mL | 112075804 |
| Lysing Matrix A | 50 ea | 116910050 | 5 ea | 116910005 |
| Lysis Buffer R | 60 mL | 116541052 | 6 mL | 116541002 |
| Wash Buffer R | 12 mL | 116543051 | 1.2 mL | 116541003 |
| Nuclease-free water | 10 mL | 116541054 | 1.5 mL | 116541004 |
| DNase I | 1 vial | 116541055 | 1 vial | 116541055 |
| DNase I Buffer | 5 mL | 116541056 | 0.5 mL | 116541006 |
| Column R with collection tube | 50 ea | 116541057 | 5 ea | 116541007 |
| Quick-start protocol | 1 ea | - | 1 ea | - |
| Instruction Manual | Available www.mpbio.com | | | |
| MSDS & CoA | Available www.mpbio.com | | | |

2.2 User Supplied Materials

- FastPrep Instrument - FastPrep-24™ 5G (Cat. No.116005500)
- Vortex mixer with adapter
- Microcentrifuge capable of at least 14,000 g
- Absolute ethanol (100 mL for preparing Wash Buffer R and 750 µL per prep for sample preparation)
- Nuclease-free 2 mL microcentrifuge tubes
- Nuclease-free 1.5 mL microcentrifuge tubes
- Single-channel pipettors (1 µL-1000 µL)
- Nuclease-free, aerosol-preventive tips
- Biohazard disposal containers
- Microcentrifuge tube rack
- Personal Protective Equipment

3. Storage and Kit Stability

DNase I should be stored at 2-8°C. All other components and reagents of the SPINeasy RNA Kit for Tissue (With Lysing Matrix) can be stored at room temperature (15-25°C) until the expiration date printed on the kit label. For extended storage or storage in dry condition (humidity < 40%), store the columns at 2-8°C to maintain performance.

4. Important Consideration Before Use

- ☐ Add 100 mL (10 mL for sample kit) of absolute ethanol to Wash Buffer R and mark the bottle.
- ☐ Prepare DNase I solution according to instructions in Section 6 (No.1).
- ☐ Prepare one 1.5 mL microcentrifuge tube per prep for elution of purified RNA.
- ☐ Prepare one 1.5 mL microcentrifuge tube for DNase I-buffer mixture.
- ☐ For some samples, lysis may be performed by vortexing the sample in a vial of Lysing Matrix A at the maximum speed if a FastPrep Instrument is unavailable.
- ☐ Centrifugation speed stated in the manual will be a guideline; use the maximum speed available if 14,000 g is not feasible.

5. Safety Precautions

Lysis Buffer R contains a component that can be harmful if swallowed and may cause irritation when in contact with skin and eyes. To prevent accidental ingestion, do not eat, drink, or smoke when using this product. Wear personal protective equipment (gloves, lab coat and eye protection) to prevent contact with the skin or mucous membranes. Consult the Material Safety Data Sheet at www.mpbio.com for additional details.

6. Protocol

1. Preparation of DNase I Solution

- Spin down briefly the vial of lyophilized **DNase I** provided and resuspend with **500 µL Nuclease-free water**. Mix well to dissolve. Store **DNase I** solution at -20 °C in aliquots and avoid repeated freeze-thawing.
- **Note:** Do not prepare **DNase I** solution in **DNase I Buffer**.

2. Homogenize

- Animal/ plant tissue: Weigh **10 - 30 mg** of animal tissue or **50 - 100 mg** of plant tissue. Cut tissue into small pieces, transfer to a vial of **Lysing Matrix A** and add **1 mL Lysis Buffer R**.
- Cell culture: Resuspend cell pellet (1×10^6 cells recommended) in **1 mL Lysis Buffer R** and transfer to a vial of **Lysing Matrix A**.
- Homogenize in a FastPrep Instrument for **15 sec @ 4.0 m/s** speed setting.
- Homogenization by FastPrep is recommended, especially for plant tissues. If a FastPrep Instrument is not available, lysis may be performed by vortexing samples in **Lysing Matrix A** for **3 - 5 min @ maximum speed**.
- Centrifuge for **10 min @ 14,000 g**.
- **Note:** Centrifuge at the maximum speed for all steps if **14,000 g** is not feasible.
- Carefully pipette out **750 µL** of the supernatant to a clean, nuclease-free **2 mL microcentrifuge tube**.

3. Bind

- Add **750 µL absolute ethanol** and mix well by pipetting up and down.
- Transfer **750 µL** of the mixture to a **Column R** with collection tube.
- Centrifuge for **1 min @ 14,000 g**. Discard flow through and reuse collection tube.
- Repeat to load the remaining mixture.

4. Wash & DNase I treatment

- Add **500 µL Wash Buffer R** to the column.
- Centrifuge for **1 min @ 14,000 g**. Discard flow through and reuse collection tube.
- **DNase I digestion:**

- In a clean 1.5 mL microcentrifuge tube, add **5 µL DNase I solution** to **75 µL DNase I Buffer** per prep. Mix well and add **80 µL** to the center of the column membrane.
- Incubate at room temperature for 15 mins.

5. Wash

- Add **500 µL Wash Buffer R** to the column.
- Centrifuge for **1 min @ 14,000 g**. Discard flow through and reuse collection tube.
- Add **500 µL Wash Buffer R** to the column.
- Centrifuge for **1 min @ 14,000 g**. Discard flow through and reuse collection tube.
- Centrifuge for additional **1 min @ 14,000 g** to dry column.

6. Elution

- Remove collection tube and place column onto a clean 1.5 mL microcentrifuge tube.
- Add **100 µL Nuclease-free water** to the center of the membrane. Incubate for **1 min @ room temperature**. For samples with low RNA content, reducing the elution volume to **50 µL** may increase the concentration of eluted RNA.
- Centrifuge for **1 - 2 min @ 8,000 g** to elute RNA.
- Eluted RNA will be collected in the microcentrifuge tube. For the best results, proceed to perform downstream applications immediately and keep RNA chilled on ice while working, to prevent degradation. Store remaining RNA at **-80 °C** in aliquots and avoid repeated freeze-thawing.
-

7. Data

The following results are examples of RNA extracted from various animal and plant tissues using SPINeasy RNA Kit for Tissue. Extracted RNA are of high quality and integrity, as indicated by A_{260}/A_{280} and A_{260}/A_{230} ratios of around 2, as well as 28S and 18S bands clearly resolved in agarose gel electrophoresis. RNA extracted with this kit can be amplified by RT-PCR.

Table 1: Quality and quantity of RNA extracted from various tissues using SPINeasy RNA Kit for Tissue.

| Sample | Starting Amount (mg) | Extraction Results | | |
|------------------------------------|-------------------------|--------------------------------|---------------|---------------|
| | | Concentration (ng/ μ L) | $A_{260/280}$ | $A_{260/230}$ |
| Rabbit Liver | 21.4 | 626.46 | 2.08 | 2.19 |
| Rabbit Spleen | 14.0 | 310.53 | 2.04 | 2.26 |
| Rabbit Kidney | 22.3 | 237.88 | 2.06 | 2.21 |
| Rabbit Lung | 21.4 | 159.24 | 2.01 | 2.26 |
| Rabbit Heart | 22.5 | 54.56 | 2.01 | 2.82 |
| Ginkgo biloba leaf | 51.3 | 132.93 | 2.09 | 2.21 |
| Prunus davidiana leaf (wild peach) | 52.7 | 382.08 | 2.07 | 1.94 |
| Cherry leaf | 48.5 | 294.06 | 2.06 | 1.91 |
| Begonia leaf | 52.0 | 368.66 | 2.06 | 1.94 |
| Tomato leaf | 52.2 | 663.82 | 2.10 | 2.28 |
| Peanut leaf | 70.2 | 898.92 | 2.10 | 2.31 |
| Potato leaf | 52.4 | 491.20 | 2.14 | 2.30 |




Figure 1: RNA extracted from various samples using SPINeasy RNA Kit for Tissue, analyzed using agarose gel electrophoresis. M: DNA marker; Lane 1: Liver; Lane 2: Spleen; Lane 3: Kidney; Lane 4: Lung; Lane 5: Heart; Lane 6: Ginkgo biloba leaf; Lane 7: Prunus davidiana leaf; Lane 8: Cherry leaf; Lane 9: Begonia leaf; Lane 10: Tomato leaf; Lane 11: Peanut leaf; Lane 12: Potato leaf.




Figure 2: RT-PCR amplification of RNA extracted from various samples using SPINeasy RNA Kit for Tissue. M: DNA marker; Lane 1: Liver; Lane 2: Spleen; Lane 3: Kidney; Lane 4: Heart; Lane 5: Lung; Lane 6: Negative control (rabbit β -actin); Lane 7: Ginkgo biloba leaf; Lane 8: Prunus davidiana leaf; Lane 9: Cherry leaf; Lane 10: Begonia leaf; Lane 11: Tomato leaf; Lane 12: Peanut leaf; Lane 13: Potato leaf; Lane 14: Negative control (plant 18S).

8. Troubleshooting

This guide may be useful in solving any problems that may arise. For further assistance, please contact our technical support team at apac-techsupport@mpbio.com

| Problem | Possible Cause | Recommendation |
|-------------------|----------------------------------|---|
| Low/No RNA Yield | Inefficient extraction | Ensure the extraction was performed according to kit's manual instructions. |
| | Sample with low RNA content | (i) Increase amount of starting material; (ii) Elute in a smaller volume (50 µL). |
| | Insufficient lysis or over-lysis | When using vortex instead of FastPrep for sample homogenization, lysis condition may be optimized by testing reduced or extended vortexing duration. For tough samples such as plant tissues, lysis by vortexing could be inefficient and result in low RNA yields. For such samples, homogenization by FastPrep is highly recommended. |
| | Poor sample quality | For best results, freshly prepared samples should be used. |
| | RNase contamination | Work with nuclease-free tubes and pipette tips. Handle samples and perform all steps with clean gloves. Decontaminate work surfaces with RNase Erase® (Cat. No. 112440204). |
| Smeared RNA bands | Poor sample quality | For best results, freshly prepared samples should be used. |
| | Sample over-lysis | Reduce FastPrep speed and/or duration. |
| | RNA degradation | Work with freshly purified RNA and keep RNA chilled on ice after elution. RNA should be stored at -80 °C, freeze thawing should be avoided. |
| DNA contamination | Inefficient DNase I digestion. | Perform on-column DNase I digestion according to step 4 of the RNA extraction protocol. Make sure that DNase I solution is prepared according to kit's manual instructions. Once dissolved, DNase I should be stored at -20 °C. |

9. Product Use Limitation & Warranty

The products presented in this instruction manual are for research or manufacturing use only. They are not to be used as drugs or medical devices to diagnose, cure, mitigate, treat, or prevent diseases in humans or animals, either as part of an accepted course of therapy or in experimental clinical investigation. These products are not to be used as food, food additives or general household items. Purchase of MP Biomedicals products does not grant rights to reproduce, modify, or repackage the products or any derivative thereof to third parties. MP Biomedicals makes no warranty of any kind, expressed or implied, including merchantability or fitness for any particular purpose, except that the products sold will meet our specifications at the time of delivery.

Buyer's exclusive remedy and the sole liability of MP Biomedicals hereunder shall be limited to, at our discretion, no replacement or compensation, product credits, refund of the purchase price of, or the replacement of materials that do not meet our specification. By acceptance of the product, Buyer indemnifies and holds MP Biomedicals harmless against, and assumes all liability for, the consequence of its use or misuse by the Buyer, its employees, or others, including, but not limited to, the cost of handling. Said refund or replacement is conditioned on Buyer notifying within thirty (30) days of receipt of product. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

FastDNA®, FastRNA®, FastPrep®, QBiogene®, and BIO 101® Systems are trademarks of MP Biomedicals, LLC.

Australia

Tel: +61 2.8824.2100
Tel: +61 1800.249.998
Email: custserv.au@mpbio.com

Austria & Germany

Tel: 0800.426.67.337
Tel: 00800.7777.9999
Email: custserv.de@mpbio.com

Belgium

Tel: 00800.7777.9999
Email: custserv.be@mpbio.com

Canada

Tel: +1 800.854.0530
Email: custserv.ca@mpbio.com

China

Tel: +86 400.150.0680
Email: custserv.cn@mpbio.com

Europe

Tel: +33 3.88.67.54.25
Tel: +33 00800.7777.9999
Email: custserv.eur@mpbio.com

France

Tel: +33 3.88.67.54.25
Email: custserv.fr@mpbio.com

India

Tel: +91 22.27636921/22/25
Email: custserv.in@mpbio.com

Italy

Tel: 00800.7777.9999
Email: custserv.it@mpbio.com

Japan

Tel: +81 3.6667.0730
Email: custserv.jp@mpbio.com

Latin America

Tel: +1 800.854.0530
Tel: +1 440.337.1200
Email: custserv.la@mpbio.com

New Zealand

Tel: +64 9.912.2460
Email: custserv.nz@mpbio.com

North America

Tel: +1 800.854.0530
Tel: +1 440.337.1200
Email: custserv.na@mpbio.com

Poland

Tel: 00800.7777.9999
Email: custserv.po@mpbio.com

Russia

Tel: +7 495 604.13.44
Email: custserv.rs@mpbio.com

Serbia

Tel: +381 11.242.1972
Email: custserv.se@mpbio.com

Singapore/ APAC

Tel: +65 6775.0008
Tel: +65 6394.7675
Email: custserv.ap@mpbio.com

South Korea

Tel: +82 2.425.5991
Email: custserv.kr@mpbio.com

Switzerland

Tel: 00800.7777.9999
Email: custserv.ch@mpbio.com

The Netherlands

Tel: 00800.7777.9999
Email: custserv.nl@mpbio.com

United Kingdom

Tel: 0800.282.474
Email: custserv.uk@mpbio.com

www.mpbio.com

