

SPINeasy™ DNA/RNA Kit for Feces

Cat. No.: 116555050 (50 PREPS)/116555000 (5 PREPS)



Scan QR code for more information from instruction manual


Quick-Start Protocol

Revision Aug 2023

This protocol is designed to extract microbial DNA and RNA from the same feces sample.

Notes before starting:

- Phenol:Chloroform:Isoamyl Alcohol (25:24:1, pH 6.5-8.0) is required for this protocol.
- The protocol requires Ethanol (96–100%) that needs to be supplied by user. Do not use denatured alcohol.
- Buffer FDR2 needs to be stored at 2–8°C upon reception.
- For faster processing, it is advised to pre-position the plasticwares and columns as shown below.



Column and sample preparation

1. Transfer the **DNA Column** (capless) into a **DNA Collection tube** (provided). Add **200 µL Equilibration Buffer** to both **DNA** (capless) and **RNA (blue) Column** membranes to ensure their performances.

Wait at least **1 min** and centrifuge for **30 sec @ maximum speed**. Transfer the Columns into new Collection tubes (provided). The treated Columns can be stored at 2–8°C for a few days if required.

2. Weight up to **250 mg** of feces and add it into **Lysing Matrix YB** tube.
3. Add **100 µL** Phenol:Chloroform:Isoamyl Alcohol and **650 µL Buffer FDR1** to the Lysing Matrix YB tube.
Homogenize using vortex at **2,500-2,700 rpm for 20 min** or **Fastprep® 35 sec, 5m/sec**. Centrifuge for **2 min @ maximum speed**.

Remove inhibitor

4. Transfer the supernatant (~400-600 µL) into a **2 mL Tube** (provided).

Optional: The **Lysing Matrix YB** tube may be re-centrifuged after transfer of the supernatant, to recover any additional lysate using a tip.

Add **200 µL Buffer FDR2** to precipitate contaminants, vortex for **1 sec** and centrifuge for **2 min @ maximum speed**.

5. DNA and RNA binding

Transfer the supernatant into a 2 mL Tube (provided). Add 1 volume **Buffer FDR3 (1:1)** and vortex for **1 sec**. For example, for 700 μL supernatant, add 700 μL Buffer FDR3.

- The treated DNA Column (capless) shall be loaded in 2 steps as follow:

- a) Pour **$\sim 750 \mu\text{L}$** of the lysate into the DNA Column (capless). Centrifuge for **10 sec @ 15,000 g**. The 2 mL Tube including the remaining lysate mixed with Buffer FDR3 can be centrifuged together to recover any mixture found on the tubing lid and wall.

- Transfer the DNA Column (capless) into a new DNA Collection tube (provided) and close the lid of the tube containing the first flow-through (it will be used for RNA isolation).

- b) Load the remaining lysate to the DNA Column (capless) and centrifuge for **10 sec @ 15,000 g**.
- Transfer the DNA Column (capless) into a new DNA Collection tube (provided) and combine the second flow-through to the first flow-through for RNA purification.

RNA purification

Add **250 μL Ethanol** (not provided) into the lysate flow-through and mix immediately by pipetting up and down or vortex for **1 sec**.

Pour **$\sim 825 \mu\text{L}$** of the obtained supernatant into the RNA Column (**blue**), centrifuge for **10 sec @ 15,000 g** and discard the flow-through.

Optional: The DNA Collection tube including the lysate-FDR3-Ethanol mix may be centrifuged together to recover any mixture found on the tubing lid and wall.

Repeat until all the lysate has been processed through the RNA Column (**blue**). Transfer the RNA Column (**blue**) to a new **RNA Collection tube** (provided).

DNA purification

If RNA and DNA is required, close the lid of the DNA Column (capless) to prevent it from drying out (until step 5 is performed). The DNA Column (capless) can also be stored at 2-8°C for extended storage.



6. **1st wash.** Add **600 μL Buffer FDR4/DNA** to the center of the DNA Column and **600 μL Buffer FDR4/RNA** to the center of the RNA Column. Centrifuge for **10 sec @ 15,000 g**, discard the flow-through and place the column back into the same Collection tube.
7. **2nd wash.** Add **600 μL Buffer FDR5** to the center of the column and centrifuge for **30 sec @ 15,000 g**, discard the flow-through and place the column back into the same Collection tube.
8. **Dry.** Transfer the column into new DNA and/or RNA Collection tube (provided) and centrifuge for **1 min @ maximum speed**.

9. Transfer the column into new **DNA and/or RNA Elution tube** (provided). Add **100 μL RNase Free Water** directly into the column membrane, wait for **1 min** and centrifuge for **1 min @ maximum speed**. The RNA/DNA samples are now ready for downstream applications.

Note: Eluting with 100 μL RNase Free Water will maximize nucleic acid yield. For more concentrated sample, a minimum of 50 μL RNase Free Water can be used.


SPINeasy™ DNA/RNA Kit for Feces

Cat. No.:11655050 (50 PREPS)/11655000 (5 PREPS)



Revision Aug 2023

Flow-Chart





MP Biomedicals


MP Biomedicals nucleic acid extraction kits are designed for simple, efficient, and rapid purification of DNA and RNA from various types of samples. Our wide range of instruments and reagent kits provide you a one-stop solution for your sample preparation works.




FastPrep® Homogenizer
Powerful, Flexible, Intuitive



Lysing Matrix
Complete, Quantitative lysing method



MagBeads FastDNA / FastRNA Kits & SPINeasy™ Extraction Kits
Highly purified DNA / RNA / Protein



Automated Nucleic Acid Purification (aNAP) Instruments
MPure-96™ Process up to 96 samples
MPure-32™ Process up to 32 samples



MP BIOMEDICALS
APAC - +65 6775 0088 | custserv.ap@mpbio.com
EUROPE - 00800 777 9999 | custserv.eu@mpbio.com
AMERICAS - 800 854 9391 | custserv.us@mpbio.com
Learn more at www.mpbio.com

