

# SPINeasy™ DNA Kit for Microbiome

Cat. No.: 116553050 (50 PREPS) / 116553000 (5 PREPS)



## Quick-Start Protocol

Revision Oct 2023

The kit allow efficient and unbiased lysis of microbes including gram positive/negative bacteria, fungi, protozoans, and viruses found in various samples including soil, swabs, body fluids, milk etc.

### Notes before starting:

- If **Buffer MB1** has precipitated, heat at **37 °C** until precipitate dissolves.
- Centrifugation speed stated in the manual will be a guideline; use the maximum speed available if **15,000 g** is not feasible.
- Vortex the samples at **2,500-3,000 rpm** for **10 min** if a **FastPrep®** instrument is unavailable. Secure samples on the vortex through an adaptor to ensure homogenization.
- For fast processing, the protocol is compatible with vacuum manifold.

Scan QR code for more information from instruction manual



Column Preparation

### Optional: Prepare column to ensure its performance

Pipette **200 µL Equilibration Buffer** into **Column MB**. Incubate for **1 min** at room temperature and centrifuge for **30 sec @14,000 g**. **Discard flow-through and reuse collection tube.**

Lyse & Remove inhibitor

### DNA isolation protocol:

1. Add **10<sup>7</sup>~10<sup>10</sup> cfu bacterial sample** to **Lysing Matrix E tube**.  
**Note:** Refer to next page for recommended amounts and summary of sample preparation.
2. Add **800 µL Buffer MB1** and **25-40 µL RNase A Solution**, invert the tubes several times to mix the lysing matrix and buffer. Homogenize in a **FastPrep®** instrument for **45 sec** twice at speed setting of **5 m/s** with a **5 minutes** interval. Centrifuge for **2 min @ ≥15,000 g**.
3. Transfer all supernatant into 2 mL centrifuge tube (self-provided). Add **300 µL Buffer MB2**, invert and mix 5 times. Centrifuge for **2 min @ ≥15,000 g**.
4. Transfer the supernatant (~750 µL) to a clean 2 mL centrifuge tube (self-provided). Add **750 µL Buffer MB3**, invert and mix twice.

Binding

### Microcentrifuge Method

5. Transfer **~750 µL** mixture to **Column MB**.
6. Centrifuge for **1 min @ ≥15,000 g** and discard the flow-through. Repeat the process until all the lysate has passed through.
7. Add **500 µL Buffer MB4** to the center of the column, centrifuge for **1 min @ 15,000 g**, discard the flow-through and place the column back into the same 2 mL collection tube.
8. Add **700 µL Buffer MB5** to the center of the column and centrifuge for **1 min @ 15,000 g** and discard the flow-through.

### Vacuum Manifold Method

5. Insert the **Column MB** into the vacuum manifold's luer connectors. Load **~750 µL** mixture into the **Column MB** and apply vacuum.
6. Repeat until all the mixture has been loaded. Switch off the vacuum source to avoid membrane over drying.
7. Add **500 µL Buffer MB4** to the center of the column and apply vacuum. Switch off the vacuum source.
8. Add **700 µL Buffer MB5** by running the pipette tip along the wall of the column and apply vacuum. Switch off the vacuum source.

Wash

9. Without addition of any liquid, centrifuge for **2 min @ 15,000 g** to dry the column.
10. Discard collection tube and place the **Column MB** into a new **1.5 mL Collection Tube** (provided).
11. Add **45-50 µL Buffer MB6** to the middle of the membrane column **slowly**. Centrifuge for **1 min @ 15,000 g**. Reload the elute or add **45-50 µL fresh Buffer MB6** into the column Centrifuge for **1 min** at maximal speed.

**Note:** For maximal yield, the elution volume can be expanded to **200 µL**.

Elute

# SPINeasy™ DNA Kit for Microbiome


Cat. No.: 116553050 (50 PREPS) / 116553000 (5 PREPS)



## Quick-Start Protocol


Revision Oct 2023

### Flow-Chart



## Summary of sample preparation in Lysing Matrix E

Sample Type	Optimal amount of sample usage	Pretreatment
Bacteria	50-100 mg (wet weight) or up to $10^{10}$ bacterial and $10^8$ yeast cells	Resuspend cell pellet in 800 $\mu$ L of Buffer MB1 and transfer to Lysing Matrix E.
Sputum	<300 $\mu$ L	Sputum added to Lysing Matrix E for lysis; for large volume sputum pretreatment, refer to manual.
Bronchoalveolar	1-10 mL	After centrifugation for 15 min @ 3,000 g (at 4 °C), recover the precipitate for lysis step.
Urine	1-20 mL	After centrifugation for 15 min @ 14,000 g, recover the precipitate for lysis step.
Milk	1-10 mL	Same as urine
Vinasse	200 mg	Add 800 $\mu$ L Buffer MB1 directly for lysis.
Soil	100-500 mg	Same as vinasse



### Pre-treatment of special samples

1. For difficult samples such as spores, increase lysis speed to 6-7 m/s; for easy samples, the lysis time can be reduced to improve efficiency.
2. Using the Buffer MB1, this kit can be used to extract DNA from oral/ nasal swabs directly and swabbing solution.
3. Refer to manual for specific sample processing methods.



MP Biomedicals

MP Biomedicals nucleic acid extraction kits are designed for simple, efficient, and rapid purification of DNA and RNA from various types of samples. Our wide range of instruments and reagent kits provide you a one-stop solution for your sample preparation works.




**FastPrep® Homogenizer**

Powerful, Flexible, Intuitive



**Lysing Matrix**

Complete, Quantitative lysing method



**MagBeads FastDNA / FastRNA Kits & SPINeasy™ Extraction Kits**

Highly purified DNA / RNA / Protein



**Automated Nucleic Acid Purification (aNAP) Instruments**

MPure-96™ Process up to 96 samples

MPure-32™ Process up to 32 samples



MP BIOMEDICALS

APAC - +65 6775 0088 | [custserv.ap@mpbio.com](mailto:custserv.ap@mpbio.com)  
EUROPE - 0800 777 9999 | [custserv.eu@mpbio.com](mailto:custserv.eu@mpbio.com)  
AMERICAS - 800 854 9393 | [custserv.us@mpbio.com](mailto:custserv.us@mpbio.com)  
Learn more at [www.mpbio.com](http://www.mpbio.com)

