

SPINeasy™ DNA Kit for Tissue (With Lysing Matrix)



Cat. No.: 116558050 (50 PREPS) / 116558000 (5 PREPS)

Quick-Start Protocol


Revision Nov 2023



Scan QR code for more information
from instruction manual

Notes before starting:

- Add 12 mL (1.2 mL for sample kit) of absolute ethanol into **Buffer TD3** and mark the bottle.
- Add 50 mL (5 mL for sample kit) of absolute ethanol into **Buffer TD4** and mark the bottle.
- This kit requires the use of a centrifuge capable of generating at least 14,000 g to obtain optimal results. Use the maximum speed available if 14,000 g is not feasible.
- If FastPrep-24™ 5G (Cat. No.116005500) is not available, the use of a vortex capable of achieving 2,500 rpm is required.
- This Kit requires the use of a ThermoMixer capable of simultaneous shaking at 1,000 rpm and heating the samples to 56 °C.
- This kit can also be used with a vacuum manifold for the bind and wash step. Please refer to the instruction manual for more details.



Column preparation

Optional: Column preparation:

Note: Column preparation is recommended when higher DNA yield is desired or when column performance is reduced after long-term storage.

1. Pipette 200 µL Equilibration Buffer into **Column S** (assembled with **Collection tube**). Incubate for 1 min at room temperature and centrifuge the column for 30 sec @ 14,000 g.
2. Keep the columns aside for later use (The treated Columns S can be stored at 2-8 °C for up to 7 days, if required).

Sample preparation

DNA isolation protocol:

1. Weigh tissue (up to 10 mg for spleen tissue, up to 30 mg for other tissues) and place in a **Lysing Matrix M** tube.
2. Add 200 µL **Buffer TD1**, 20 µL **Proteinase K** and 4 µL **RNase A** into the tissue sample tube, vortex for 5 sec to mix well. Briefly spin down the mixture.
3. Homogenize using FastPrep® for 5 sec @ 4 m/s or vortex for 5 min @ 2,500 rpm. Briefly spin down the lysate.
4. Incubate in a ThermoMixer at 1,000 rpm for 10 min at 56 °C. Briefly spin down the lysate.

Bind

5. Add 500 µL **Buffer TD2** into the lysate. Mix thoroughly by pipetting up and down for 10 times or vortex for 10 sec. Briefly spin down the mixture.
6. Assemble **Column S** onto a clean **Collection tube**.
7. Load all the mixture (~700 µL) into **Column S**. Centrifuge for 30 sec @ 14,000 g. Discard flow through and place the column back into the same **Collection tube**.

Wash

8. Add 500 µL **Buffer TD3** onto the center of the column, centrifuge for 30 sec @ 14,000 g. Discard flow through and place the column back into the same **Collection tube**.
9. Add 500 µL **Buffer TD4** onto the center of the column, centrifuge for 30 sec @ 14,000 g. Discard flow through and place the column back into the same **Collection tube** (**Repeat this step once**).
10. Transfer the column to a **new** **Collection tube** and spin for 2 min @ maximum speed.

Elute

11. Transfer the column to a **Elution tube**. Add 50-100 µL **Buffer TD5** onto the center of the column, wait for 2 min and centrifuge for 2 min @ 14,000 g. Purified DNA is now ready for downstream applications.

Optional: Perform a second elution step with further 50-100 µL **Buffer TD5** will increase yields by up to 20%.

SPINeasy™ DNA Kit for Tissue (With Lysing Matrix)



Cat. No.: 116558050 (50 PREPS) / 116558000 (5 PREPS)

Revision Nov 2023

Flow-Chart

Column
preparation

Sample preparation

Bind

Wash

Elute

Optional:
Incubate 1 min at
room temperature
🌀 14,000 g, 30 sec



Equilibration Buffer
200 µL



Weigh tissue samples and
add them into a tube of
Lysing Matrix M

4 m/s
5 sec



Or

Buffer TD1 200 µL
Proteinase K 20 µL
RNase A 4 µL
Mix well

🌀 Quick spin

2,500 rpm
5 min

🌀 Quick spin

56 °C, 1,000 rpm, 10 min

🌀 Quick spin

Buffer TD2 500 µL
Mix well

🌀 Quick spin

Load all the mixture
into Column S with
Collection tube

🌀 14,000 g, 30 sec

Buffer TD3 500 µL

🌀 14,000 g, 30 sec

Buffer TD4 500 µL

🌀 14,000 g, 30 sec

Transfer the column to a
new Collection tube

🌀 maximum speed, 2 min
Column drying

Transfer the column
to Elution tube

Buffer TD5 50-100 µL
Incubate at RT for 2 min

🌀 14,000 g, 2 min

Highly purified
genomic DNA



MP BIOMEDICALS

APAC +65 6775 0068 | custserv.ap@mpbio.com
EUROPE +30 800 777 9999 | custserv.eu@mpbio.com
AMERICAS +00 858 9293 | custserv.us@mpbio.com
Learn more at www.mpbio.com

