

Rapid SARS-CoV-2 Antigen Test Card

CE 0123 Certified Self-Test for Private Use

The Rapid SARS-CoV-2 Antigen Test Card is an in-vitro diagnostic test for private use. It is used to detect SARS-CoV-2 antigens from anterior nasal swabs.

- ▶ **EASY TO USE:** Easy sampling using anterior nasal swab specimens
- ▶ **RAPID:** Results in 15 minutes
- ▶ **CONVENIENT:** Works with a nasal swab of 2.5 cm
- ▶ **EXCELLENT PERFORMANCE CHARACTERISTICS:** Specificity 99.20%, Sensitivity 96.77%, Accuracy 98.72%
- ▶ **PACKAGE:** Includes all required consumable materials
- ▶ **NO CROSS-CONTAMINATION** with other coronaviruses like MERS, 229E, NL63, OC43 or influenza A or B



Pack Size	1 Test/Box (Cat. No.: 07AG6001BS) 5 Tests/Box (Cat. No.: 07AG6005BS) 20 Tests/Box (Cat. No.: 07AG6020BS) includes all necessary materials
Application	Identification of SARS-CoV-2 infection
Specimen	Nasal / anterior nasal
Duration	Results after 15 minutes
Interpretation	Results are easy to read
Storage Temperature	4 - 30 °C



ORDER NOW!



MP Biomedicals Germany GmbH

✉ laientest@mpbio.com

☎ 0049.5651.921111

Learn more at www.mpbio.com



SARS-CoV-2AG-SELF-FL2102EN