



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

ImmunOTM

Catalog Number: 421651

Immunodiffusion Plates, Uncalibrated, unfilled

Description: Plates are polystyrene with close-fitting covers to minimize gel dehydration. Plate bottoms are recessed 7 mm to prevent scratching on the bench surfaces. The plate is hydrophilic for easy gel pouring and better adherence to the plate. The working surface of the plate is 95 mm by 45 mm.

Use: MP's immunodiffusion plates are intended for use in double immunodiffusion and single radial immunodiffusion.

Calibration: These plates are uncalibrated.

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.