

MP Biomedicals, LLC

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone: 440/337-1200  
Toll Free: 800/854-0530  
Fax: 440/337-1180  
mailto: [biotech@mpbio.com](mailto:biotech@mpbio.com)  
web: <http://www.mpbio.com>

## TECHNICAL INFORMATION

Catalog Number: 151944, 191427, 194726, 194844

### Potassium Chloride

**Molecular Formula:** KCl

**Molecular Weight:** 74.55

**CAS #** 7447-40-7

**Physical Description:** White crystalline powder

**Solubility:** Soluble in water (1 g/3 ml - clear, colorless solution), boiling water (1 g/1.8 ml), glycerol (1 g/14 ml), or ethanol (1 g/250 ml); insoluble in ether or acetone; hydrochloric acid, sodium or magnesium chlorides diminish its solubility in water.<sup>1</sup>

**Description:** A salt commonly used in the preparation of phosphate buffered saline and for the extraction and solubilization of proteins.<sup>3,5</sup> A good potassium source.<sup>1,4</sup> Also used in the hypotonic KCl-propidium iodide protocol for univariate analysis of metaphase chromosomes<sup>2</sup>; dissociation of eukaryotic ribosomes into subunits by high concentrations of potassium (0.3-1.0 M)<sup>6</sup>; preparation of mitochondria, provides an ionic osmotic support for those tissues that assume a gelatinous consistency upon homogenization.<sup>7</sup>

### Availability:

Catalog Number	Description	Size
151944	Potassium Chloride, purity approximately 99%	250 g 500 g 1 kg 5 kg
194726	Potassium Chloride, cell culture reagent, purity approximately 99%	250 g 500 g 1 kg 5 kg
194844	Potassium Chloride, molecular biology reagent, purity approximately 99%	500 g 1 kg 5 kg
191427	Potassium Chloride, ACS Reagent Grade	100 g 500 g 1 kg 5 kg

### References:

- *Merck Index, 12th Ed.*, No. 7783.
- Cram, L.S., et al., *Meths. Cell Biol.*, v. **33**, 369 (1990).
- Garcia-Hernandez, M., et al., *J. Biol. Chem.*, v. **269**, 20744-20749 (1994).
- Henshaw, E.C. and Panniers, R., *Meths. Enzymol.*, v. **101**, 616 (1983).
- Huvos, P. and Cox, R.A., *Biochim. Biophys. Acta*, v. **383**, 421-426 (1975).
- Labowitz, A.M., *Meths. Enzymol.*, v. **59**, 421 (1979).
- Nedergaard, J. and Cannon, B., *Meths. Enzymol.*, v. **55**, 3 (1979).