



MP Biomedicals, LLC

29525 Fountain Parkway  
Solon, Ohio 44139

Telephone: 440/337-1200  
Toll Free: 800/854-0530  
Fax: 440/337-1180  
mailto: [biotech@mpbio.com](mailto:biotech@mpbio.com)  
web: <http://www.mpbio.com>

## TECHNICAL INFORMATION

### ENZYME SYSTEMS PRODUCTS

a division of



### Technical Data Sheet

#### General Instructions for Use of Chloromethyl Ketone (CMK) Products

##### Storage:

In solid form, the material may be stored at room temperature in a desiccator. In DMSO or DMF solution, store at -20° C.

##### Solubility:

Highly soluble in DMSO and DMF (1-100 mM)

##### Use:

Dilute from stock into enzyme/buffer solution after enzyme pre-incubation/activation. Allow to act, as necessary (optimal incubation time to be determined by researcher). Ten to twenty minutes of inhibitor incubation is suggested.

##### References:

– Segal, D.M. et al (1971). *Biochem* **10**: 3728