



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

Catalog Number: 1004149

EY Tellurite Enrichment

To be used with Baird-Parker Agar Base (Catalog # 10-015) for isolating *S. aureus*.

Composition: Consists of 30% egg yolk suspended with 0.15% potassium tellurite.

Use: Fifty ml of EY Tellurite Enrichment is added aseptically to 950 ml of prepared and sterilized Baird-Parker Agar Base which has been cooled to 50°C.

Caution: EY Tellurite Enrichment, because of its nature, cannot be sterilized by usual procedures. The treatment given this product cannot guarantee the destruction of heat resistant spores of bacteria or yeasts. This product has passed our sterility test and use and performance tests. However, it is recommended that the usual sterility test performed on the completed medium be read before it is used.