



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337 , 1200
Toll Free: 800/854 , 0530
Fax: 440/337 , 1180
mailto:biotech@mpbio.com
web: www.mpbio.com

TECHNICAL INFORMATION

Catalog Number: 64140

Bovine IgG

Content : Lyophilized powder at a pH of between 4 and 5

Procedures : IgG is prepared from bovine serum using current chromatographic techniques. A quantitative protein analysis as well as gel - electrophoresis and immunoassay techniques have been applied to determine quality of the product.

Reconstitution : Reconstitute in 5.0 ml of water or 0.15 M NaCl for a 2% solution

Purity : Characterization by gel-electrophoresis

Note : This product does not contain any preservative such as sodium azide, thimerosal or proclin. Neither does the product contain any surfactants