



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

Catalog Number: 151121

Ferritin

Physical Description: Red to Yellow Solution

Source: Human Liver

Form: Suspension in 0.15M NaCl containing 0.1% sodium azide

Protein Content: ≥ 2.25 mg/ml

Purity: Not less than 95% as determined by PAGE

Solubility: Soluble in distilled water or dilute buffer

Warning: Human source material. Treat as potentially infectious. Each donor unit used in the preparation of this product has been tested by and FDA approved method and found non-reactive for the presence of I IBsAg (an antibody to HIV). Because no known test method can offer complete assurance that hepatitis B virus, HIV, or other infectious agents are absent, all human blood-based products should be handled in accordance with good laboratory practices using appropriate precautions.*

*Centers of Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories", 1988.