



MP Biomedicals, LLC

29525 Fountain Parkway
Solon, Ohio 44139

Telephone: 440/337-1200
Toll Free: 800/854-0530
Fax: 440/337-1180
mailto: biotech@mpbio.com
web: <http://www.mpbio.com>

TECHNICAL INFORMATION

cappel

Catalog Number: 55876

Washed, Preserved Sheep Red Blood Cells (10% Suspension)

DESCRIPTION: Product is 20ml of a 10% suspension of sheep red blood cells collected into Alsever's solution (Dextrose, Sodium Chloride, Sodium Citrate), washed and preserved with 1% Penicillin/Streptomycin.

USES: This product is suitable for immunoassays (rosetting), hemolytic and hemagglutination assays and for assays of heme inhibition.

SPECIFICATIONS: Product reacts in an agglutination assay with anti-sheep red blood cells. The hematocrit is within the acceptable range.

RELATED PRODUCT ORDERING INFORMATION:

Catalog #DescriptionSize

55800 RABBIT ANTI-SHEEP RED BLOOD CELLS2ML

55806 IGG RABBIT ANTI-SHEEP RED BLOOD CELLS5ML

55808 IGM RABBIT ANTI-SHEEP RED BLOOD CELLS5ML

Note: This product may contain a preservative such as sodium azide, thimerosal or proclin. Please see lot specific chemical credential for preservative information.